Health o meter[®]
P r o f e s s i o n a l

Model 450KL Physician Beam Scale

ISO 9001:2015
CERTIFIED

ISO13485:2016
CERTIFIED

User Instructions

PELSTAR, LLC 9500 West 55th St. McCook, IL 60525-7110 USA www.homscales.com

Health o meter[®] Professional 450KL

© Pelstar, LLC 2025

weigh**easier**[®]

Thank you for purchasing this Health o meter® Professional product. Please read this manual carefully, and keep it for easy reference or training.

TABLE OF CONTENTS

Cautions and Warnings	3
Specifications / Certifications / Service	4
Assembly Instructions	5
Operating Instructions	7
Height Measurement Instructions.....	9
Counterweight Instructions.....	11
Transportation	15
Troubleshooting	16
Maintenance and Cleaning.....	18
Warranty.....	19

Note: This scale has been factory calibrated, and does not require calibration prior to use. There are no requirements for calibration in the field; users should adhere to their institution’s calibration policies.

Please register your scale for warranty coverage at:

www.homscales.com

For User Instructions updates and revisions please go to:

www.homscales.com

CAUTIONS AND WARNINGS

INTENDED USE

This Health o meter® Professional platform scale is intended to be used in a professional medical environment by trained medical staff. This product was designed to weigh patients who are safely positioned and standing on the platform. The intended use of the built-in height rod is to measure patient height. Do not modify the product or use it for anything other than its intended purpose.

To prevent patient/caregiver injury or damage to your scale, please follow the instructions in this user manual very carefully.

- Do not exceed the specified weight limit for this scale.
- Do not transport the scale with a patient or object on the scale.
- Do not wheel or pull the scale down stairs, doing so may damage internal parts.
- To prevent injury, as well as scale damage during assembly, exercise caution when assembling the scale pillar.
- Assemble and operate the scale per the enclosed user instructions.
- For accurate weighing, this scale must be placed on a flat, stable surface.
- For accurate weighing, verify before each use the proper operation according to the procedure described in this manual.
- Do not use in the presence of flammable or explosive materials.
- If the scale becomes damaged, it should not be operated until properly serviced.

Patient/caregiver safety

- This scale is designed for static weighing of patients only. No scale should be used for patient transfer.
- Patient should wear socks or other form of lightweight foot covering during weighing event.
- To prevent patient injury, the patient must be attended to throughout the entire weighing event. Caregivers should ensure that the patient is stable and provide support as needed when getting onto and off of the scale.

In no event whatsoever shall Pelstar, LLC be liable for damages or injuries arising from or connected with the assembly, use, or misuse of its products.

SPECIFICATIONS / CERTIFICATIONS / SERVICE

General

The Model 450KL physician beam scale employs a proven and rugged mechanical design to measure patient weight and height. Each precision scale is designed to provide accurate, reliable, and repeatable weight measurements. In addition, each scale is designed with a rotating poise bar, and includes counterweights that enable the user to make the weighing process simple, fast, and convenient.

Weight is displayed in pounds (lb) or in kilograms (kg). Height is displayed in inches (in) or in centimeters (cm).

450KL Scale Specifications

Capacity:	350 lb (160 kg)	
Capacity with included counterweights:	500 lb (200 kg)	
Environmental:	Operating temperatures: 50°F to 104°F (10°C to 40°C) Storage temperatures: 32°F to 122°F (0°C to 50°C) Maximum Humidity: 85% RH	
Physical Dimensions:	Scale (Overall Dimensions) Width: 10-1/2" (26.7 cm) Depth: 20-1/4" (51.4 cm) Height: 57-1/2" (146.1 cm) Weight: 28 lb (12.7 kg)	Platform (Standing Area) Width: 10-1/2" (26.7cm) Depth: 14" (35.6 cm) Height: 3-1/4" (8.3 cm)

Certification Descriptions

 	<p>Pelstar, LLC has been officially certified as the manufacturer of Health o meter® Professional medical devices. The Pelstar, LLC quality assurance system covers the development, production, sales and service of medical scales and measuring systems.</p>
--	---

Customer Service Information

The 450KL physician beam scale is shipped disassembled in one carton. Carefully inspect the carton for shipping damage before unpacking. If damage is found, contact Customer Service.

For further information or telephone support, please contact Customer Service at:

Telephone: 1 (800) 815-6615

Email: homprocs@hom scales.com

ASSEMBLY INSTRUCTIONS

Before Assembly

Each 450KL physician beam scale is shipped disassembled in one carton. Carefully inspect the carton for shipping damage before unpacking. If damage is found, contact your shipper or a Health o meter® Professional representative immediately. Claims must be filed with the shipper as soon as possible after receipt of the package. The following information details what you will find inside the main carton as you unpack the parts for assembly.

To prevent scratching any components and avoid damaging the scale parts when unpacking, do not use a box cutter, knife, scissors, or any sharp object to open the protective inner packaging. Carefully remove each assembly from the carton and unwrap the packing materials. Set the carton aside for storage.

Parts List

Carton

- (1) Scale platform
- (1) Pillar with height rod attached
- (4) Pillar nuts
- (1) User instructions
- (1) Wrench
- (1) Counterweight set (Part Number 55070)
 - 150 lb counterweight (Part Number 55060)
 - 40 kg counterweight (Part Number 55065)


Tools Required

Wrench (included)


ASSEMBLY INSTRUCTIONS (CONT.)

Caution: It is recommended that the assembly be performed by two people. Use of a chair is also recommended for the assembly process.

1. Carefully remove the scale platform and pillar from the carton and place on a flat, level, and dry surface. Locate the steel rod at the base of the pillar. Cut the tie wrap that holds the steel rod / finger pull loop to the pillar.


2. Obtain a sturdy chair and place the top of the pillar with the height rod side facing down on the seat of the chair, as shown. Place the platform/base perpendicular to the platform with the front of the base resting on the floor and the top of the platform/base facing the top of the pillar. Align the four studs on the bottom of the pillar with the four holes in the top of the platform/base. Insert the studs into the platform/base.


Note: If a chair is not available it is recommended that one person holds the platform/base while another person handles the pillar.


3. Access the underside of the scale base. Obtain the four pillar nuts in the carton and place onto the studs on the underside of the scale base. Using the provided wrench, tighten the pillar nuts to secure the pillar to the base.

4. Grasp the finger pull loop on the end of the steel rod and pull towards the lever extension. Push the lever extension towards the finger pull loop. Hook the steel rod to the hole at the end of the lever extension. Release the finger pull loop and ensure that the steel rod and lever extension are hooked together. Carefully position the scale to its upright position. The scale is now assembled. Refer to the next page to zero balance the scale before first use.

OPERATING INSTRUCTIONS

Zero Balance Adjustment

Before first use, the scale needs to be zero balanced by following the instructions below. The scale should also be checked and balanced periodically, particularly when moved.

Note: A flat head screwdriver is required.

1. Ensure that the scale is located on a hard, flat, and level surface. Step on and off scale several times in order to align all internal parts.
2. Place both poise weights on the "0" value on both bars. The beam pointer should float gently up and down, and not touch the top or bottom of the trig-loop.
3. If the beam pointer is not centered in the trig-loop opening, using a flat head screwdriver, rotate the balance ball screw to correct the beam pointer. Turn the screw about $\frac{1}{4}$ turn forwards or backwards and let poise bar settle. Based on which direction the poise bar moved turn $\frac{1}{4}$ forwards or backwards again and let poise bar settle until the pointer is centered in the trig loop.

Weighing Instructions

1. Set the upper and lower poise weights to indicate the patient's approximate weight.
2. Have patient step onto the scale. Ensure they are positioned directly in the center of the platform. Continue to adjust the upper and lower poise weights until the pointer is centered.



Note: It is not necessary for you to wait for the pointer to stop moving in order to read the correct weight. As long as the pointer is moving equally within the trig-loop opening, the patient's weight can be determined.

3. The patient weight can be recorded in either pounds or kilograms. Rotate the upper and lower poise bars to measure in either pounds or kilograms. For easy identification, pounds (lb) have silver graduations on a black background, located on one side of the poise bar; kilograms (kg) have black graduations on a silver background, located on the other side of the poise bar. Ensure that poise bars are locked after being rotated into position.
4. Observe the location of the indicator arrow built into each poise weight. Read the value on the poise bar, closest to the indicator arrow.
5. Add the indicated reading of the upper poise value to the indicated reading of the lower poise value. The total amount is the patient's weight. See the example on the next page.

Note: To transport your scale, acquire the optional wheel bracket (Item # 55000). Please refer to the "Transportation" section on page 9 for assembly instructions.

Note: To increase your scale's capacity, utilize the included counterweight set. Please refer to the "Counterweight Instructions" section on page 10.

OPERATING INSTRUCTIONS (CONT.)

HEIGHT MEASUREMENT INSTRUCTIONS

WARNING: Use caution when sliding height rod up and down. The rod should move smoothly, do not move the rod in a jerking motion. Fold headpiece down when not in use.

Note: A patient needs to be a minimum of 11 7/8" (30 cm) tall for height to be measured.

1. With the headpiece in the folded down position, slide the headpiece up or down to the estimated height of the patient. Then raise the headpiece into its horizontal position.

NOTE: FOLD HEADPIECE DOWN BEFORE EXTENDING HEIGHTROD

Push tab inward to lower height rod below 46 3/8"

To lower height rod below 46 3/8", push the tab at the top of the height rod and gently lower the headpiece.

2. Ensure that the patient's posture is upright, and that the patient's head is straight and level. Slide the headpiece down until it rests on the patient's head.



3. Determine the patient's height by reading the value on the measuring strip nearest to the line on the indicator, as shown on the next page.

HEIGHT MEASUREMENT INSTRUCTIONS (CONT.)

Height Rod Below 46 3/8"

Height Rod Extended Above 46 3/8"

Indicator

Inches read: 41" + 3/8" = 41 3/8"


Cm read: 105 cm + 0.2 cm = 105.2 cm

Inches read: 69" + 3/8" = 69 3/8"

Cm read: 176 cm + 0.2 cm = 176.2 cm

Feet and Inches read: 3' + 5 3/8" = 3 ft. 5 3/8"

Feet and Inches read: 5' + 9 3/8" = 5 ft. 9 3/8"

COUNTERWEIGHT INSTRUCTIONS

In order to increase the scale's capacity to a maximum value of 500 lb / 200 kg, the counterweight set must be used.

To employ the counterweight set, refer to the instructions below. Contact technical support with any questions or concerns regarding the use of the counterweight set.


Note: Use counterweight 55060 for pounds and counterweight 55065 for kilograms.

Counterweight Instructions for Weighing in Pounds (lb)

Note: The 150 lb (55060) counterweight is required when weighing a patient who is over 350 lbs.

Counterweight Instructions for Weighing a Patient 351 - 388 (lb)

1. When a patient weighs between 351 - 388 lb, place the lower poise weight at 250 lb.
2. Place the 150 lb (55060) counterweight on the 3rd notch from the left, between the 50 lb and 100 lb value.
3. Adjust the upper poise weight to match the patient weight. The poise bar will center in the trig-loop.
4. Read and record the weight. To read the patient's weight, add the weight value indicated on the upper and lower poise bars, as well as an additional 88 lb. See page 13 for a counterweight conversion chart.

Counterweight Instructions for Weighing a Patient 389 - 424 (lb)

1. When a patient weighs between 389 - 424 lb, place the lower poise weight at 300 lb.
2. Place the 150 lb (55060) counterweight on the 3rd notch from the left, between the 50 lb and 100 lb value.
3. Adjust the upper poise weight to match the patient weight. The poise bar will center in the trig-loop.
4. Read and record the weight. To read the patient's weight, add the weight value indicated on the upper and lower poise bars, as well as an additional 88 lb. See page 13 for a counterweight conversion chart.

COUNTERWEIGHT INSTRUCTIONS (CONT.)

Counterweight Instructions for Weighing a Patient 425 - 450 (lb)

1. When a patient weighs between 425 - 450 lb, place the lower poise weight at 250 lb.
2. Place the 150 lb (55060) counterweight on the 5th notch from the left, between the 150 lb and 200 lb value.
3. Adjust the upper poise weight to match the patient weight. The poise bar will center in the trig-loop.
4. Read and record the weight. To read the patient's weight, add the weight value indicated on the upper and lower poise bars, as well as an additional 150 lb. See page 13 for a counterweight conversion chart.

Counterweight Instructions for Weighing a Patient 451 - 500 (lb)

1. When a patient weighs between 451 - 500 lb, place the lower poise weight at 300 lb.
2. Place the 150 lb (55060) counterweight on the 5th notch from the left, between the 150 lb and 200 lb value.
3. Adjust the upper poise weight to match the patient weight. The poise bar will center in the trig-loop.
4. Read and record the weight. To read the patient's weight, add the weight value indicated on the upper and lower poise bars, as well as an additional 150 lb. See page 13 for a counterweight conversion chart.

COUNTERWEIGHT INSTRUCTIONS (CONT.)

Counterweight Instructions for Weighing in Kilograms (kg)


Note: The 40 kg (55065) counterweight is required when weighing a patient who is over 160 kg.

Counterweight Instructions for Weighing a Patient 161 - 180 (kg)

1. When a patient weighs between 161 - 180 kg, place the lower poise weight at 120 kg.
2. Place the 40 kg (55065) counterweight on the 5th notch from the left, between 60 kg and 80 kg value on the lower poise bar.
3. Adjust the upper poise weight to match the patient weight. The poise bar will center in the trig-loop.
4. Read and record the weight. To read the patient's weight, add the weight value indicated on the upper and lower poise bars, as well as an additional 40 kg. See page 13 for a counterweight conversion chart.

Counterweight Instructions for Weighing a Patient 181 - 200 (kg)

1. When a patient weighs between 181 - 200 kg, place the lower poise weight at 140 kg.
2. Place the 40 kg (55065) counterweight on the 5th notch from the left, between 60 kg and 80 kg value on the lower poise bar.
3. Adjust the upper poise weight to match the patient weight. The poise bar will center in the trig-loop.
4. Read and record the weight. To read the patient's weight, add the weight value indicated on the upper and lower poise bars, as well as an additional 40 kg. See page 13 for a counterweight conversion chart.

COUNTERWEIGHT INSTRUCTIONS (CONT.)

Counterweight Conversion Chart





Counterweight Set 55070															
Counterweight 55060 - Pounds					Counterweight 55065 - Kilograms										
Combine the weight shown on the upper and lower poise then add 88lbs.					Combine the weight of the upper and lower poise then add 150lbs.					Combine the weight shown on the upper and lower poise then add 40Kg.					
Counterweight between 50 & 100 Lower poise weight at 250		Counterweight between 50 & 100 Lower poise weight at 300			Counterweight between 150 and 200 Lower poise weight at 250		Counterweight between 150 and 200 Lower poise weight at 300			Counterweight between 60 & 80 Lower poise weight at 120		Counterweight between 60 & 80 Lower poise weight at 140			
If the scale reads	The patient's weight is	If the scale reads	The patient's weight is	If the scale reads	The patient's weight is	If the scale reads	The patient's weight is	If the scale reads	The patient's weight is	If the scale reads	The patient's weight is	If the scale reads	The patient's weight is		
262	350	301	389	275	425	301	451	121	161	141	181	121	161	141	181
263	351	302	390	276	426	302	452	122	162	142	182	122	162	142	182
264	352	303	391	277	427	303	453	123	163	143	183	123	163	143	183
265	353	304	392	278	428	304	454	124	164	144	184	124	164	144	184
266	354	305	393	279	429	305	455	125	165	145	185	125	165	145	185
267	355	306	394	280	430	306	456	126	166	146	186	126	166	146	186
268	356	307	395	281	431	307	457	127	167	147	187	127	167	147	187
269	357	308	396	282	432	308	458	128	168	148	188	128	168	148	188
270	358	309	397	283	433	309	459	129	169	149	189	129	169	149	189
271	359	310	398	284	434	310	460	130	170	150	190	130	170	150	190
272	360	311	399	285	435	311	461	131	171	151	191	131	171	151	191
273	361	312	400	286	436	312	462	132	172	152	192	132	172	152	192
274	362	313	401	287	437	313	463	133	173	153	193	133	173	153	193
275	363	314	402	288	438	314	464	134	174	154	194	134	174	154	194
276	364	315	403	289	439	315	465	135	175	155	195	135	175	155	195
277	365	316	404	290	440	316	466	136	176	156	196	136	176	156	196
278	366	317	405	291	441	317	467	137	177	157	197	137	177	157	197
279	367	318	406	292	442	318	468	138	178	158	198	138	178	158	198
280	368	319	407	293	443	319	469	139	179	159	199	139	179	159	199
281	369	320	408	294	444	320	470	140	180	160	200	140	180	160	200
282	370	321	409	295	445	321	471								
283	371	322	410	296	446	322	472								
284	372	323	411	297	447	323	473								
285	373	324	412	298	448	324	474								
286	374	325	413	299	449	325	475								
287	375	326	414	300	450	326	476								
288	376	327	415			327	477								
289	377	328	416			328	478								
290	378	329	417			329	479								
291	379	330	418			330	480								
292	380	331	419			331	481								
293	381	332	420			332	482								
294	382	333	421			333	483								
295	383	334	422			334	484								
296	384	335	423			335	485								
297	385	336	424			336	486								
298	386					337	487								
299	387					338	488								
300	388					339	489								
						340	490								
						341	491								
						342	492								
						343	493								
						344	494								
						345	495								
						346	496								
						347	497								
						348	498								
						349	499								
						350	500								

TRANSPORTATION

A wheel bracket is included if purchased as model 450KLWH. If purchased as model 450KL an optional wheel bracket can be purchased separately (Item #55000). Follow these steps to install the wheel bracket:

Tools Required (not included)

Adjustable wrench or 10 mm wrench

<p>1. Remove the mounting screws that attach the clamping plate to the wheel bracket.</p>	 <p>Wheel bracket</p> <p>Clamping plate</p> <p>Mounting screws</p>
<p>2. Gain access to the underside of the scale base by carefully tipping the scale forward until the pillar is resting on the floor. Position the clamping plate inside the base by aligning the holes in the clamping plate and base.</p>	 <p>Mount wheel bracket here</p> <p>Close-up view of wheel mounting area</p> <p>Base</p>
<p>3. Position the wheel bracket onto the base by aligning the holes in the wheel bracket with the holes in the clamping plate.</p>	
<p>4. Insert the mounting screws through the wheel bracket and into the clamping plate. Using a wrench, tighten the mounting screws to secure the wheel bracket to the base. Carefully tip the scale to its upright position.</p>	 <p>Mounting screws</p>

To transport the scale, tip the pillar backwards until the scale is resting on the wheel bracket. Carefully move the scale to its new location.

TROUBLESHOOTING

Troubleshooting

Before contacting service personnel, refer to the following instructions to check and to correct any failures for incorrect weight measurements.

For the most accurate weight measurement, use the scale on a hard level surface and stand in the center of the platform. If the error seems excessive, check the following:


Symptom	Possible Cause	Corrective Action
Questionable weight, or the scale does not center equally in the trig loop when the poise weight is "0" value	Bar does not "0" balance	The balance ball must be balanced such that the pointer comes to rest in the center of the trig-loop when both poise weights are set at "0" value. Ensure that the upper and lower poise weights are firmly against the shoulder of the beam. Adjust the balance ball screw by turning the screw at the left end of the bar until balanced.
Bar does not move freely	Bar sticking or adjustment needed	Check if the pointer is touching the side of the trig loop in its range of travel. Realignment of the beam may be necessary if it is touching. Ensure that the upper and lower poise weights are firmly against the shoulder of the beam. Adjust the balance ball screw by turning the screw at the left end of the bar.
Platform rocks excessively or touches base at any corner	Base out of alignment	See instructions on next page.
Beam does not move at all during weighing process	Scale may have improper steel rod connection or poise weight may have been set at a higher weight than patient's actual weight.	Reset the weights to a lower weight. Ensure the steel rod is connected properly, refer to "Assembly Instructions" on page 5.
Scale not accurately measuring weight		Contact Customer Service at 1 (800)-638-3722.

For further assistance please contact technical support toll free at 1 (800) 638-3722.

TROUBLESHOOTING (cont.)

Base Alignment Instructions


1. Check the space between the platform and the base. The space should be even on all four sides of the base. To test if the space is even gently press down on one of the toe-end corners and note the depressed space. Press down on the other toe-end corner and compare the space. Do the same for all four corners to detect which corner is higher or lower than others.


2. To align the space between the platform and base the hex screws on the underside of the base need to be adjusted. To access the underside of the base obtain a sturdy chair and place the top of the pillar with the height rod side facing down on the seat of the chair, as shown.

3. On the underside of the scale base remove the two black caps from the "heel end".

4. With the caps removed you can access the alignment screws that allows the "heel end" of the base to move up or down. Use a $7/64$ " or $1/8$ " allen wrench to adjust the screws. As you adjust the screws, watch the space between the base and platform move to ensure the space is even. When the space looks even, test the platform by pressing each corner one time. The space on all four corners should be the same as the "toe end". If all are even replace the black caps and return the scale to its upright position.

MAINTENANCE AND CLEANING

The following provide instructions for maintenance and cleaning for the 450KL. Maintenance operations other than those described in this manual should be performed by qualified service personnel.

Maintenance

Before first use and after periods of non-use, check the scale for proper operation and function. If the scale does not operate correctly, refer to qualified service personnel.

1. Check overall appearance of the total scale for any obvious damage, wear and tear.
2. Do not store near chemicals, heaters or radiators.
3. Do not drop, kick, or jump on the scale.
4. This scale is a sensitive weighing device. Items such as clothing or towels placed or dropped on the scale may cause it to operate incorrectly. In order to prevent this, do not store or leave anything on the scale.
5. If your scale fails to work properly, do not attempt to dismantle it. This will void your warranty. See warranty information instructions on having your scale serviced.
6. All bearings and working parts are permanently lubricated. **DO NOT USE OIL OR LUBRICANTS.**

Cleaning and Disinfecting

Proper care and cleaning is essential to ensure a long life of accurate and effective weighing.

1. Health o meter® Professional recommends using one of the following solutions on a soft cloth or disposable wipe:
 - mild soap and water solution
 - 70% isopropyl alcohol
 - solution with 1-5% hydrogen peroxide concentrationAfter cleaning/disinfecting, wipe with a cloth dampened with water and then with a clean dry cloth. To avoid residue build-up or damp surfaces, ensure all scale parts are completely dry after cleaning.
2. Never use rough or abrasive materials to clean the scale, as these will damage the scale's finish.
3. Do not submerge the scale in water or any other liquid.
4. Do not pour or spray fluids directly on the scale.

Replacement Parts

Replacement parts may be available if a piece of the scale needs to be replaced. Contact Health o meter Professional Scales Customer Service at 1-800-815-6615 to inquire about the availability of these replacement parts.

PART NUMBER	DESCRIPTION	UNIT OF MEASURE
400HARDWARE	ALLEN WRENCH AND HEX SCREWS	1 SET
55000	WHEELS	1 EA
55070	COUNTERWEIGHT SET (LB AND KG)	2 EA
402KLROD	HEIGHT ROD	1 EA
400BASE	BASE	1 EA
400COVER	PLATFORM COVER	1 EA

WARRANTY

10-Year Limited Warranty

What does the Warranty Cover?

This Health o meter® Professional scale is warranted from date of purchase against defects of materials or in workmanship for a period of ten (10) years. If product fails to function properly, return the product, freight prepaid and properly packed to Pelstar, LLC (see "To Get Warranty Service", below, for instructions). If the manufacturer determines that a defect of material or in workmanship exists, the customer's sole remedy will be replacement of the scale at no charge. Replacement will be made with a new or remanufactured product or component. If the product is no longer available, replacement may be made with a similar product of equal or greater value. All replaced parts are covered only for the original warranty period.

Who is Covered?

The original purchaser of the product must have proof of purchase to receive warranty service. Please save your invoice or receipt. Pelstar dealers or retail stores selling Pelstar products do not have the right to alter, or modify or any way change the terms and conditions of this warranty.

What is Excluded?

Your warranty does not cover normal wear of parts or damage resulting from any of the following: negligent use or misuse of the product, use on improper voltage or current, use contrary to the operating instructions, abuse including tampering, damage in transit, or unauthorized repair or alternations. Further, the warranty does not cover natural disasters, such as fire, flood, hurricanes and tornadoes. This warranty gives you specific legal rights, and you may also have other rights that vary from country to country, state to state, province to province or jurisdiction to jurisdiction.

To get Warranty Service make sure you keep your sales receipt or document showing proof of purchase. Call (+1) 800-638-3722 or (+1) 708-377-0600 to receive a return authorization (RA) number, which must be included on the return label. Attach your proof of purchase to your defective product along with your name, address, daytime telephone number and description of the problem. Carefully package the product and send with shipping and insurance prepaid to:

Pelstar, LLC
Attention R/A# _____
Return Department
9500 West 55th Street
McCook, IL 60525

PELSTAR, LLC

9500 West 55th St. McCook, IL 60525-7110 USA

1-800-638-3722 or 1-708-377-0600

PLEASE REGISTER YOUR SCALE FOR WARRANTY COVERAGE AT:

www.homscales.com

Health o meter® is a registered trademark of Sunbeam Products, Inc. used under license. Health o meter® Professional products are manufactured, designed, and owned by Pelstar, LLC. We reserve the right to improve, enhance, or modify Health o meter® Professional product features or specifications without notice.

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PELSTAR®, LLC
© Pelstar®, LLC 2025