

Powered by

VALHALLA
SCIENTIFIC

Body Composition Scale

ADVANTAGES

✔ **15 Second Measurement Time**

On screen display delivers rapid access to weight, body fat %, body water %, lean mass, Basal Metabolic Rate, & Body Mass Index. A full illustrated printout can be generated in a minute

✔ **Large Platform**

The spacious platform allows patients to find a comfortable and secure stance.

✔ **Few Pretest Requirements**

Analysis can be performed at any time of day and under almost any circumstances, resulting in the greatest flexibility of any body composition analyzer.

✔ **High Capacity**

The scale features an industry leading 1000 lb capacity, that provides the flexibility and accuracy to weigh any patient.

✔ **Patient Report Templates***

The G6 Series scales offer users 4 graphic printout formats to use with their patients: Adult/Bariatric, Adult Spanish, Pediatric and Fitness.

* User must connect scale to a printer to provide personalized graphic printout. Printer is optional and sold separately, not available through Health o meter® Professional scales.

 **ScaleSurance Certified**

Extended Warranty Program that offers an extra 2 years of coverage

Powered by

VALHALLA
SCIENTIFIC

Body Composition Scale

Model	Shipping Dims. (L x W x H) / Cube	Weight (LB)	UPC Code
BCS-G6	46" x 19" x 10" / 10	37 lb	854896003702

Specifications

- Capacity: 1000 lb / 453.5 Kg
- Resolution: 0.1 lb / 0.1 kg
- Platform Size: 14" (w) x 16" (d) x 1.875" (h)
- Display: 6-window display (1/2" digits) includes weight, % body fat, % body water, % muscle mass, BMI, REE
- Power Source: 110-220V power adapter (included)
- Optional printer, sold separately, not available through Health o meter® Professional scales
- 1 Year Limited Warranty - Eligible for ScaleSurance Extended Warranty

Bio-Impedance Specifications

- Electrode Method: Tetrapolar System
- Frequency: 50kHz
- Test Current: 500µA
- Conversion Rate: 750msec



BCS-G6-adult Body Comp Scale	The whole-body, adult only scale is ideal for practices that exclusively see patients 18 and over	854896003702
BCS-G6-limbs Body Comp Scale	The adult segmented version provides patients with body composition information about their torso as well as specific information about the composition of their arms and legs	854896003719
BCS-G6-DUO Body Comp Scale	The whole-body adult/pediatric version is a good choice for providers who see adults and children	854896003733
BCS-G61-upper Body Comp Scale	The adult upper-body-only configuration provides a partial body composition analysis for patients who do not want to remove their shoes and socks	854896003726
IPO-ADU	Spanish Language Patient Report Template (5-page report, 25 reports/package)	854896003740
IPO-PED	Spanish Language Patient Report Template (5-page report, 25 reports/package)	854896003757
IPO-ATH	Spanish Language Patient Report Template (5-page report, 25 reports/package)	854896003764
IPO-ADU-SP	Spanish Language Patient Report Template (5-page report, 25 reports/package)	854896003948
BCS-BROCHURE	Waiting Room Patient Brochure	