σ О S S Φ O Ω



Advantages

330 Ib Capacity Floor scale

The high capacity of this floor scale provides more versatile use

Durable Construction

Solid construction and durable, non-skid mat can withstand frequent use

Easy-to-read Dial

Large, easy-to-read dial is conveniently located at the top of the scale

Small-Step Platform

Low platform height requires only minimal step up, offering patient a more comfortable and stable weighing experience



See reverse side for additional specifications





Specifications

- Capacity: 330 lb
- Graduation: 1 lb
- Platform Size (w x d x h): 10 ¹/₂" x 9 ³/₆" x 2" / 266 mm x 238 mm x 51 mm
- Product Footprint (w x d x h): 11 ¾" x 9 ¾" x 2" / 289 mm x 248 mm x 57 mm
- Optional Carrying Case Available (order # 64771)
- 1 Year Limited Warranty



Carrying Case for 175LB Item# 64771



Health o meter® Professional Scales 9500 West 55th Street USA McCook, IL 60525-7110 Phone: 1-800-815-6615 Fax: 1-800-654-7330

(ഷ)

B

0

S

S

Ф

 \bigcirc

175LB_F&B_20191021