



## Advantages

- ✓ **Moveable Arm & Foot Rest**  
The flip-arm and foot rest allows for easier patient access and patient comfort
- ✓ **Motion Detection**  
Scale has the ability to weigh accurately despite movement of the subject on the scale
- ✓ **Value-Added Features**  
Operating functions including LB/KG Conversion, LB/KG Lock Out, Zero Out, Tare, Pre-Tare, Auto Off, Reweigh and more!



Model	Shipping Dims. (L x W x H) / Cube	Weight (LB)	UPC Code
2595KL	32.625" x 24.875" x 41.875" / 19.67	100.25 lb	892076002626



**Specifications**

- Capacity: 600 lb / 272 kg
- Resolution: 0.2 lb / 0.1 kg
- Seat Size: 18 1/4"(w) x 15 3/4" (d) x / 464 mm (w) x 400 mm (d)
- Display: 1 1/8" LCD
- Wheels: 2
- Connectivity: via optional Pelstar® wireless technology
- Power Source: 110-240V Adapter (included), 6 C-cell batteries (not included)
- Functions: LB / KG Conversion, LB / KG Lock Out, Zero Out, Tare, Auto Off, Variable auto-off Time, Enabling unit view with unit of measure locked, Hold/Release, Reweigh, Pre-Tare, Recall, Everlock®
- 2 Year Limited Warranty

Model	Description
2595KL	Digital Chair Scale - pounds & kilograms
2595KG	Digital Chair Scale - kilograms only
ADPT30	Optional adapter for 2595 series of scales
PLUGSET30	Plug adapters for use outside North America

Scale also includes connectivity capabilities. To learn about connectivity options, visit [homsscales.com/ConnectivitySolutions](http://homsscales.com/ConnectivitySolutions)