

386S01 (POUNDS)



386KGS01 (KILOGRAMS)



Advantages

- ✓ **Sanitary Plastic Weighing Surface**
The easy-to-clean plastic tray is ideal for sterile environments
- ✓ **Built-In Measuring Tape**
Easy-to-read measuring tape allows user to measure both weight and height simultaneously
- ✓ **Portable**
Convenient carrying handle for easy transport

386KGS01 (KILOGRAMS)



386S01 (POUNDS)



Model	Master Pack	Shipping Dims. (L x W x H) / Cube	Weight (LB)	UPC Code
386S01	1	16.25" x 21.5" x 6.88" / 1.39	11 lb	892076002305
386KGS01	1	16.25" x 21.5" x 6.88" / 1.39	11 lb	892076002312

Specifications

- Capacity: 50 lb (386S01)
- Capacity: 25 kg (386KGS01)
- Graduation: 1/4 lb (386S01)
- Graduation: 100 g (386KGS01)
- Tray Size: 19-3/8" (w) x 12-1/8" (d) x 2-1/8" (h)
- Easy to clean plastic tray
- Measuring Tape: 0" - 18" / 0 cm - 46 cm
- 1 Year Limited Warranty