

## THE 522KL-HR

# DIGITAL PEDIATRIC TRAY SCALE AND MECHANICAL HEIGHT ROD

The 522KL-HR pediatric scale was specifically created for nurseries, delivery rooms, and pediatric practices, making it an ideal choice for newborns, infants, and toddlers. With a maximum capacity of 50 pounds, this scale offers providers and patients a range of benefits, including safety, comfort, and flexibility. The easy to read height rod makes measuring infant height fast and easy. The scale ensures a secure platform and keeps patients calm for accurate readings.

The 522KL-HR display is positioned 17" above the scale base, ensuring an unobstructed easy-to-read weight measurement. The user-friendly digital interface includes features such as Zero, Tare, Pre-Tare, Hold/Release, and Auto-off, making it convenient and easy to use. Its motion detection technology ensures accurate weight measurements even when the patient moves.

The 522KL-HR provides easy integration and seamless communication with today's EMR applications with the option to incorporate wireless connectivity for interfacing with computers, EMR software, or other electronic devices.



## THERE'S STRENGTH IN OUR NUMBERS.

When you add up the numbers on our scale, it's not surprising that the 522KL is your choice for safety, comfort, and reliability.

**120+** LOCAL  
FIELD SPECIALISTS

**30-60 SEC**  
AUTO WEIGHT HOLD

**#1** TRUSTED  
BRAND

**100 YEARS**  
OF SERVICE

**Health o meter®**  
Professional

Health o meter® Professional 522KL-HR

weigh**easier**®

# TECHNICAL SPECIFICATIONS

## Capacity

50 lb / 23 kg

## Resolution

0 - 20 lb / 0.2 oz; 20 - 50 lb / 0.5 oz;

0 - 9 kg / 5g; 9 - 23 kg / 10 g

## EMR Connectivity

Via optional Pelstar®  
wireless technology

## Display

1" / 25 mm, LCD display

## Power Source

6 AA batteries included,  
100-240V adapter (ADPT31)  
not included

## Functions

LB / KG conversion

LB / KG lock out

Zero

Tare

Pre-Tare

Hold / release

Reweigh

Auto zero

Auto off

Motion detection

Audible/mute sound option

Variable auto-off time

Enabling unit view with unit of measure locked

Auto-weight hold for 30 to 60 seconds

Everlock® - permanent weight lock

## Tray Size (w x d x h)

24 1/8" x 14 1/2" x 2 5/8"

(991 mm x 368 mm x 67 mm)

## Product Footprint

(height rod extended)

39" x 21" x 23 1/8"

(800 mm x 533 mm x 587 mm)

## Product Weight

20 lb / 9 kg

## Mechanical Height Range

3" - 39" / 7 cm - 100 cm

Graduation: 1/16" / 1 mm

## Warranty

2-year limited manufacturer warranty

# ORDERING INFORMATION

MODEL	DESCRIPTION
522KL-HR	522KL with attachable mechanical height rod
522KG-HR	522KG with attachable mechanical height rod
522HR	Mechanical height rod for 522 series of scales
SS-522KL	ScaleSurance® 2 year extended warranty for 522 series*
ADPT31	Optional adapter for 522 series of scales
PLUGSET30	Plug adapters for use outside North America
BTPRINTERKIT-01	Wireless thermal printer kit

Scale also includes connectivity capabilities. To learn about connectivity options, visit [homscales.com/ConnectivitySolutions](https://homscales.com/ConnectivitySolutions)

To learn more, call **1-800-815-6615**  
or visit [homscales.com/522KL-HR](https://homscales.com/522KL-HR)

**Health o meter®**  
Professional



\*Not available in all countries.