

Health o meter®

P r o f e s s i o n a l

Large 1 ½" LCD Display



Advantages

- ✓ **1 ½" LCD Display**
The large display shows the user an easy-to-read weight measurement
- ✓ **400 lb / 180 kg Capacity Floor Scale**
The high capacity of this floor scale provides more versatile use
- ✓ **5-Second Weight Hold**
Weight measurement will stay displayed for 5 seconds after stepping off scale allowing user time to stabilize themselves or the patient



Carrying Case for 800KL
Item# 64771



Specifications

- Capacity: 400 lb / 180 kg
- Resolution: 0.1 lb / 0.05 kg
- 1-1/2" LCD display
- Functions: LB / KG Switch, Auto Zero, 5 Second Weight Hold, Auto Off, Tap-On to Start
- Platform Size (w x d x h):
12 1/4" x 12 1/4" x 1 3/4" / 311 mm x 311 mm x 44 mm
- Product Weight: 4 lbs / 1.8 kg
- (2) AAA batteries included
- Optional Carrying Case Available (order # 64771)
- 1 Year Limited Warranty