



Advantages

- ✓ **500 lb / 220 kg Capacity Floor Scale**
This high capacity floor scale is versatile for most patients
- ✓ **Large 2-1/2" LCD Display**
The large display offers easy-to-read weight measurement
- ✓ **Wide Platform**
The spacious platform allows patients to find a comfortable and secure stance to help maintain their balance
- ✓ **10-Second Weight Hold**
Weight measurement stays displayed for 10 seconds after stepping off the scale so patients can stabilize themselves
- ✓ **Low Profile Platform**
Low platform height requires only minimal step up, offering patient a more stable weighing experience
- ✓ **Non-Slip Surface**
The platform's textured mat provides traction and a stable foothold



Optional Carrying Case Available for 880KLS Order Item #64769

Model	Master Pack	Shipping Dims. (L x W x H) / Cube	Weight (LB)	UPC Code
880KL	4	18.5" x 14.4" x 6.9" / 1.06	42 lb	854896003665
880KLS	1	20" x 14" x 3"	10 lb	



MADE IN HUNGARY
LATEX FREE

Specifications

- Capacity: 500 lb / 220 kg
- Resolution: 0.2 lb / 0.1 kg
- Platform Size: 17" (w) x 13" (d) x 1-3/10" (h) / 432 mm (w) x 330 mm (d) x 33 mm (h)
- Display: 2-1/2"
- Power Source: (2) 3-volt lithium batteries (included)
- Functions: LB / KG / ST Switch, Auto Off, Tap-On to Start
- 10 Second Weight Hold
- Optional Carrying Case for 880KLS (Item # 64769)
- 1 Year Limited Warranty