

THE EBABYHR

DIGITAL BABY HEIGHT ROD

The Health o meter® Professional eBABYHR is a portable digital height rod for measuring infants. Its lightweight, and easy-to-read with a digital display that allows users to view measurements in centimeters, feet, or inches. The large head and foot paddles allow for fast precise measurement of infants while lying down.



TECHNICAL SPECIFICATIONS

Measuring Range

13 3/4" – 31 1/2" / 35 cm – 80 cm

Resolution

1/16" / 0.1 cm

Product Footprint (w x d x h)

(height rod closed)

24 3/4" x 11 1/2" x 2 3/4"

(630 mm x 290 mm x 70 mm)

Display

LCD, 5/8" / 16 mm digit height

Product Weight

1.4 lb / 0.6 kg

Power Source

2 AAA batteries (included)

Functions

Auto-Off

Hold

Warranty

2-year limited manufacturer warranty

ORDERING INFORMATION

MODEL	DESCRIPTION
eBABYHR	Digital Baby Height Rod
522EHR	eBABYHR that attaches to 522KL & 522KG
553EHR	eBABYHR that attaches to 553KL & 553KG
550EHRBRACKET	550KL mounting bracket for eBABYHR (must purchase with eBABYHR to mount on 550KL)
565EHR	eBABYHR that attaches to 565KL & 565KG
2210EHR	eBABYHR that attaches to 2210KL-AM & 2210KG-AM**
SS-eBABYHR	ScaleSurance® 2 year extended warranty for eBABYHR*

*Not available in all countries.

**Available For Use with 2210KL-AM Series 2 Sided Tray Only

Health o meter®
Professional

Health o meter® Professional eBABYHR

weigh**easier**®

To learn more, call **1-800-815-6615**
or visit **hom scales.com/eBABYHR**