

THE 2900KL-AM

ANTIMICROBIAL WALL-MOUNTED FOLD-UP WHEELCHAIR SCALE

The number-one-selling medical scale company in North America offers an antimicrobial wall-mounted fold-up wheelchair scale that supports up to 1000 pounds and provides the optimum in quality, reliability, and ease of use. The low profile scale extends off the wall 3.5" when folded up and is excellent for facilities that need to save space. The large easy-to-lift platform accommodates most wheelchairs and provides a comfortable weighing experience for patients. The dual ramps offer easy on, easy off access and the anti-slip caution tape delivers enhanced safety. Other safety features include antimicrobial protection, a raised edge, a recessed handle, and no pinch points, which provide both the clinician and patient confidence during the weighing process.

Designed to meet the needs of hospitals, clinics, acute, long-term care facilities, doctor's offices, and anywhere space is limited. The scale gives providers flexibility with a wall-mounted display that can be mounted on either side of the scale. The 2900KL-AM's high-contrast color screen provides a user-friendly interface that includes BMI, hold/ release, reweigh, recall, zero, KG/LB, tare, pretare and more. The exclusive Everlock® function allows the measuring unit (LB or KG) to be permanently locked on the scale.

The scale provides easy integration and seamless communication with today's EMR applications. The 2900KL-AM includes a USB interface for connectivity with a computer, EMR software, or other electronic devices via various Health o meter® Professional Connectivity Solutions.

THERE'S STRENGTH IN OUR NUMBERS.


When you add up the numbers, it's not surprising that the Health o meter® Professional Scales is the industry-leading name in medical scales.

120+ LOCAL FIELD SPECIALISTS

100% CONNECTIVITY

#1 TRUSTED BRAND

100 YEARS OF SERVICE



Health o meter®
Professional

Health o meter® Professional 2900KL-AM

weigh**easier**®

TECHNICAL SPECIFICATIONS

Capacity

1000 lb / 454 kg

Resolution

0.2 lb / 0.1 kg

EMR connectivity

Via USB, optional Pelstar® wireless technology

Display

Digit Size: 3/8" / 10 mm
Screen Size: 2 3/4" x 1 1/2" High-Contrast
Color TFT-LCD Screen

Power Source

100-240V adapter (ADPT30) included, 6 D batteries not included

Functions

LB/KG conversion
LB/KG lock out
BMI
Zero
Tare
Pre-tare
Hold/release
Reweigh
Recall
Auto zero
Auto off
Time/date
Variable auto off time
Audible/mute sound option
Parent/child weighing function
Everlock® - permanent unit lock

Platform Size (w x d x h)

29 1/2" x 32 1/2" x 2 3/4"
(749 mm x 819 mm x 70 mm)

Ramp Size (w x d)

6 1/4" x 34"
(159 mm x 864 mm)

Product Footprint Unfolded (w x d x h)

58" x 37" x 37"
(1473 mm x 940 mm x 940)

Product Footprint Folded (w x d x h)

58" x 3 1/2" x 37"
(1473 mm x 89 mm x 940 mm)

Product Weight

91 lb / 41 kg

Warranty

2-year limited manufacturer warranty

ORDERING INFORMATION

| MODEL | DESCRIPTION |
|-----------------|--|
| 2900KL-AM | Antimicrobial wall-mounted fold-up wheelchair scale, LB & KG |
| 2900KG-AM | Antimicrobial wall-mounted fold-up wheelchair scale, KG only |
| 2900KL-AM-BT | 2900KL-AM with built-in Pelstar® wireless technology |
| 2900KG-AM-BT | 2900KG-AM with built-in Pelstar® wireless technology |
| SS-2900KL-AM | ScaleSurance® 2 year extended warranty for 2900 series* |
| ADPT30 | Adapter included, item # for replacement or spare |
| PLUGSET30 | Plug adapters for use outside North America |
| BTPRINTERKIT-01 | Wireless thermal printer kit |

Scale also includes connectivity capabilities. To learn about connectivity options, visit homscalcs.com/ConnectivitySolutions



Easy to Lift Platform with a Recessed Handle

Health o meter®
Professional



*Not available in all countries.