

THE 2210KL-AM

ANTIMICROBIAL HIGH-RESOLUTION DIGITAL NEONATAL/PEDIATRIC TRAY SCALE

The 2210KL-AM neonatal scale was designed for use in the NICU, nursery, delivery room, and pediatric practices, where a high degree of accuracy is critical to the smallest of patients. The scale gives providers and patients safety with an Antimicrobial body, tray, handles, and keypad. The smooth-edged removable tray has a patent-pending mounting system that eliminates breakage, making it easy to clean. The built-in measuring tape allows users to measure both weight and height simultaneously. The scale will quickly lock despite the movement of the infant. The backlit display makes reading a patient's weight easier at almost any angle and light condition.

The high 1 gram resolution scale features a capacity of 20,000 grams/20 kg/45lb.


The 2210KL-AM's user-friendly digital interface includes breast milk intake, auto-weight hold, hold/release, zero, tare, pre-tare, reweigh, recall, GM/KG/LB conversion and lockout, auto-zero, and auto-off. The scale's motion detection gives it the ability to weigh accurately despite the movement of the patient. The exclusive Everlock[®] function allows the measuring units (LB, KG, or G) to be permanently locked on the scale.

2210KL4-AM Scale with 4-Sided Tray Option

The scale provides easy integration and seamless communication with today's EMR applications. The 2210KL-AM includes a USB interface for connectivity with a computer, EMR software, or other electronic devices via various Health o meter[®] Professional Connectivity Solutions. Options such as kilograms only, grams only, optional digital height rod, and built-in wireless models are available.

Health o meter[®]
Professional

Health o meter[®] Professional 2210KL-AM



THERE'S STRENGTH IN OUR NUMBERS.

When you add up the numbers on our scale, it's not surprising that the 2210KL-AM is your choice for accuracy and reliability


**1 GRAM
RESOLUTION**

**120+ LOCAL
FIELD SPECIALISTS**

**100%
CONNECTIVITY**

**#1 TRUSTED
BRAND**

**100 YEARS
OF SERVICE**



weigh**easier**[®]

TECHNICAL SPECIFICATIONS

Capacity

45 lb / 20 kg / 20,000 gm

Resolution

1 g (to 10,000 g);
5 g (10,001 to 20,000 g)

EMR Connectivity

Via USB, optional Pelstar® wireless technology

Backlit LCD Display

5 1/4" x 1 1/8" (133 mm x 28 mm)
4 levels of brightness

Power Source

6 D batteries included,
adapter ADPT30 included

Functions

LB / KG / GM conversion
LB / KG / GM lock out
Zero
Tare
Pre-tare
Reweigh
Recall
Auto zero
Auto off
Motion detection
Everlock®
Hold / release
Variable auto-off Time
Audible/mute sound option
Disable/enable backlight for display
Enabling unit view with unit of measure locked
Auto weight hold for 30 to 60 seconds
Breast milk intake

Tray Size (w x d x h)

(2210KL-AM) 24" x 15 1/4" x 2 3/4"
(610 mm x 388 mm x 71 mm)
(2210KL4-AM) 26 1/2" x 15 1/4" x 2 3/4"
(674 mm x 388 mm x 71 mm)

Product Footprint

(2210KL-AM) 24 3/4" x 17 3/4" x 6"
(630 mm x 450 mm x 153 mm)
(2210KL4-AM) 26 1/2" x 17 3/4" x 6"
(674 mm x 450 mm x 153 mm)

Product Weight

(2210KL-AM) 28 lb / 13 kg
(2210KL4-AM) 29 lb / 13.5 kg

Measuring Tape

0" - 24" / 0 cm - 60 cm

Portability

Heavy duty handles for ease of moving
2210CART-AC - optional acute care cart
2210CART - optional rolling cart

Warranty

2-year limited manufacturer warranty

ORDERING INFORMATION

MODEL	DESCRIPTION
2210KL-AM	Antimicrobial high resolution digital neonatal / pediatric tray scale, LB & KG
2210KG-AM	Antimicrobial high resolution digital neonatal / pediatric tray scale, KG only
2210KL-AM-BT	2210KL-AM with built-in Pelstar® wireless technology
2210KG-AM-BT	2210KG-AM with built-in Pelstar® wireless technology
2210KL4-AM	Antimicrobial high resolution digital neonatal / pediatric 4 sided tray scale, LB & KG
2210KG4-AM	Antimicrobial high resolution digital neonatal / pediatric 4 sided tray scale, KG only
2210KL4-AM-BT	2210KL4-AM with built-in Pelstar® wireless technology
2210KG4-AM-BT	2210KG4-AM with built-in Pelstar® wireless technology
2210CART-AC	Multi-Function Acute Care Cart for 2210KL-AM & 2210KL4-AM Series of Scales
2210CART	Rolling Scale/ Medical Device Cart
2210EHR	Digital Height Rod for the 2210KL-AM & 2210KG-AM**
SS-2210KL	ScaleSurance® 2 year extended warranty for 2210 series*
ADPT30	Adapter included, item # for replacement or spare
PLUGSET30	Plug adapters for use outside North America
BTPRINTERKIT-01	Wireless thermal printer kit

Scale also includes connectivity capabilities. To learn about connectivity options, visit homscales.com/ConnectivitySolutions

To learn more, call 1-800-815-6615
or visit homscales.com/2210KL-AM

2210EHR Optional Digital Height Rod



2210CART-AC Multi-Function Acute Care Cart for 2210KL-AM & 2210KL4-AM Series of Scales

Health o meter®
P r o f e s s i o n a l

*Not available in all countries.

**Available For Use with 2210KL-AM Series 2 Sided Tray Only