

# THE 245EHR-3001

## DIGITAL HEIGHT ROD FOR THE 3001 SERIES OF SCALES

The Health o meter® Professional digital 245EHR-3001 is an optional state-of-the-art digital height rod for the Antimicrobial 3001 series of scales. The market-leading height rod provides the most accurate height measurement in its class, with a graduation/resolution of 1 mm/0.1" and a full measuring range of 11 ¾" – 79 ½" / 30cm – 202cm and a digital range of 3ft 11.2" – 6ft 7.5" / 120 cm - 202 cm.



### Built-In Reset Feature

Proprietary technology eliminates mechanical sensors, reducing failures, increasing accuracy and reliability

The new accessory uses solid-state technology, which provides a facility confidence that a patient's height will be measured accurately, and workflow will not be interrupted. The built-in reset feature increases accuracy and reliability, and the folding headpiece limits the height rod's profile, providing safety to patients and caregivers. The optional digital height rod has quick and easy to assemble brackets.



## TECHNICAL SPECIFICATIONS

### Digital Height Range

3ft 11.2" – 6ft 7.5" / 120 cm - 202 cm  
Resolution: 1 mm / 0.1"

### Mechanical Height Range

11 ¾" – 79 ½" / 30cm – 202cm  
Graduation: 1 mm / 1/16 inch

### Product Footprint with 3001KL-AM

(height rod extended)  
22" x 33 ½" x 82 ½"  
(559 mm x 841 mm x 2096)

### Product Footprint with 3001KL-AM

(height rod in home position)  
22" x 33 ½" x 50"  
(559 mm x 841 mm x 1270 mm)

### Product Weight with 3001KL-AM

46 lb / 21 kg

### Warranty

2-year limited manufacturer warranty

### 245EHR-3001 Compatible Scales\*

3001KL-AM  
3001KG-AM  
3001KL-AM-BT  
3001KG-AM-BT

## ORDERING INFORMATION

MODEL	DESCRIPTION
245EHR-3001	Optional digital height rod for 3001 Series
SS-245EHR	ScaleSurance® 2 year extended warranty**

\*Compatible with 3001 series scales with the letter "E" or "U" in the serial number and firmware versions 02.08.00 or higher. Not available for AMX series scales. For assistance, contact Customer Service at 800-815-6615.

\*\* Not available in all countries

**Health o meter®**  
Professional

Health o meter® Professional 245EHR-3001

weigh**easier**®

To learn more, call **1-800-815-6615**  
or visit **homscales.com/245EHR-3001**