

500KL Family Comparison Guide



	500KL	500KL-BT	500KLHB	500KLAD	502KL	502KL-BT
Capacity	660 lbs / 300 kg					
Resolution	0.2 lbs / 0.1 kg					
Product Footprint	14 3/8" x 21 1/8" x 53 5/8" / 359 mm x 537 mm x 1362 mm					
Platform	13 3/4" x 16 1/2" x 2 3/8" / 349 mm x 419 mm x 60 mm					
Functions	LB/ KG Conversion, LB/ KG Lock (Temporary Lock), BMI, Zero, Tare, Pre-Tare, Hold/ Release, Auto Zero, Auto Weight Lock, Auto Off with Disable Option, Reweigh, Recall, Parent/ Child Weighing Function, Everlock® (Permanent Unit Lock)					
Mechanical Height Rod Range	76 - 214 cm/ 30 - 84.25 in	76 - 214 cm/ 30 - 84.25 in	76 - 214 cm/ 30 - 84.25 in	76 - 214 cm/ 30 - 84.25 in	30 - 202 cm/ 11.75 - 79.5 in	30 - 202 cm/ 11.75 - 79.5 in
Mechanical Height Rod Graduation	1 mm/ 1/8 in	1 mm/ 1/8 in	1 mm/ 1/8 in	1 mm/ 1/8 in	1 mm/ 1/16 in	1 mm/ 1/16 in
Digital Height Rod Range	-	-	-	-	120 - 202 cm/ 47.2 - 79.5 in	120 - 202 cm/ 47.2 - 79.5 in
Digital Height Rod Graduation	-	-	-	-	1 mm/ 1/10 in	1 mm/ 1/10 in
Wireless Connectivity	Optional	●	Optional	Optional	Optional	●
Welch Allyn Vital Signs Monitor Connectivity	Optional	●	Optional	Optional	Optional	●
Schiller Cardio, Edan/MDPro Vital Signs Monitor Connectivity	Optional	-	Optional	Optional	Optional	-
EMR Interface: Allscripts, MidMark IQ Manager, Constant Care, Charm Health, CEIBA (Direct to EMR)	Optional	●	Optional	Optional	Optional	●
License Free Software	●	●	●	●	●	●
Product Weight	21 lb/ 9.3 kg	21 lb/ 9.3 kg	31 lb/ 15 kg	21 lb/ 9.3 kg	22.4 lb/ 10.2 kg	22.4 lb/ 10.2 kg
Battery (6 AA)	●	●	●	●	●	●
Power Adapter (ADPT31)	Optional	Optional	Optional	●	Optional	Optional
EMR Capable	Optional	●	Optional	Optional	Optional	●
Transport Casters	●	●	●	●	●	●
Dent/ Rust Proof Platform	●	●	●	●	●	●
KG & LB Versions	●	●	●	●	●	●
Live Handrails (500HB)	Optional	Optional	●	Optional	Optional	Optional
Rigid Height Rod Headpiece for Better Accuracy	●	●	●	●	●	●
Meets ISO Quality Standard for Tipping Risk	●	●	●	●	●	●
No Disassembly Required for Calibration	●	●	●	●	●	●

Scale also includes connectivity capabilities. To learn about connectivity options, visit homscales.com/ConnectivitySolutions