

## THE 2210CART-AC

# MULTI-FUNCTION ACUTE CARE CART FOR 2210KL-AM SERIES OF SCALES

- 2 easy grip handles for maneuverability
- Customized plates for secure scale placement
- 2 storage shelves including a built-in 21 ¼"(w) x 13"(d) shelf and a 21 ¼"(w) x 1 ½"(d) hanging shelf
- Four large anti-static casters with locks for easy mobility
- Used to transport neonatal scale or other medical devices and supplies
- Mobile Acute Care cart designed explicitly for the Neonatal/ Pediatric scales 2210KL-AM/2210KL4-AM Series of Scales



2210CART-AC

### 2210CART-AC Compatible Scales

2210KL-AM  
2210KG-AM  
2210KL-AM-BT  
2210KG-AM-BT  
2210KL4-AM  
2210KG4-AM  
2210KL4-AM-BT  
2210KG4-AM-BT  
2210KL-AM-C  
2210KG-AM-C  
2210KL-AM-BT-C  
2210KG-AM-BT-C

## THE 2210CART

# ROLLING SCALE/ MEDICAL DEVICE CART

- The top shelf has a handle for easy grip
- 4 casters (2 casters with locks)
- Spacious storage shelf
- Lockable cabinet
- Used to transport pediatric scales or other medical devices and supplies



2210CART

### 2210CART Compatible Scales

553KL  
553KG  
549KL  
550KL  
565KL  
565KG  
565KL-C  
565KG-C  
2210KL-AM  
2210KG-AM  
2210KL-AM-BT  
2210KG-AM-BT  
2210KL4-AM  
2210KG4-AM  
2210KL4-AM-BT  
2210KG4-AM-BT  
2210KL-AM-C  
2210KG-AM-C  
2210KL-AM-BT-C  
2210KG-AM-BT-C

**Health o meter**<sup>®</sup>  
Professional

Health o meter<sup>®</sup> Professional 2210CART-AC & 2210CART

weigh**easier**<sup>®</sup>

# TECHNICAL SPECIFICATIONS

## 2210CART-AC MULTI-FUNCTION ACUTE CARE CART FOR 2210KL-AM SERIES OF SCALES

**Capacity**  
77 lb / 35 kg

**Product Footprint (w x d x h)**  
30" x 22" x 38"  
762 mm x 559 mm x 965 mm

**Storage Shelf (w x d x h)**  
21 ¼" x 13" x 3"  
540 mm x 330 mm x 76 mm

**Hanging Shelf (w x d x h)**  
21 ¼" x 1 ½"  
540 mm x 38 mm

**Wheels**  
4 large anti-static casters with locks

**Warranty**  
2-year limited manufacturer warranty

## 2210CART ROLLING SCALE / MEDICAL DEVICE CART

**Capacity**  
200 lb / 90 kg

**Product Footprint (w x d x h)**  
26 ½" x 18" x 34"  
673 mm x 457 mm x 867 mm

**Cabinet Dimensions (w x d x h)**  
23" x 16 ¾" x 17"  
584 mm x 427 mm x 432 mm

**Wheels**  
4 casters (2 with locks)

**Warranty**  
1-year limited manufacturer warranty



**2210CART-AC**  
Multi-Function Acute Care  
Cart for 2210KL-AM/  
2210KL4-AM Series of  
Scales

**2210CART**  
Rolling Scale/ Medical  
Device Cart



# ORDERING INFORMATION

MODEL	DESCRIPTION
2210CART-AC	Multi-Function Acute Care Cart for 2210KL-AM/2210KL4-AM Series of Scales
2210CART	Rolling Scale/ Medical Device Cart*

To learn more, call **1-800-815-6615**  
or visit [hom scales.com/2210CART-AC](https://hom scales.com/2210CART-AC)  
[hom scales.com/2210CART](https://hom scales.com/2210CART)



\*2210CART is a TAA compliant product.

**Health o meter®**  
Professional