

THE 565KL-C

HIGH-RESOLUTION DIGITAL BABY SCALE, CANADA ONLY

The Health o meter® Professional 565KL-C neonatal/pediatric scale is designed to meet the demands of healthcare professionals in critical settings such as the NICU, nursery, delivery room, and pediatric practices. Its exceptional accuracy is especially crucial when caring for the tiniest patients.

The gently curved cradle provides comfort for infants, while the large angled display ensures caregivers can easily read weight measurements. Adjustable feet guarantee a level weighing surface, and the built-in measuring tape allows simultaneous weight and height assessment. Enhanced motion detection ensures quick locking even with an active infant. With a high 1-gram resolution, the scale accommodates a capacity of 20 kg (44 lb).

The 565KL-C's user-friendly digital interface offers a suite of features that enhance versatility in managing neonatal/pediatric patients. Notable functions include reweigh, hold/release, tare, recall, and a 5-minute weight hold. Additionally, the scale includes zero, pre-tare, KG/LB conversion, KG/LB/GM lockout, auto-zero, and auto-off features. The exclusive Everlock® function allows the measuring units (LB or KG) to be permanently locked on the scale.

The scale provides easy integration and seamless communication with today's EMR applications. The 565KL-C offers the option to incorporate wireless connectivity for interfacing with computers, EMR software, or other electronic devices via various Health o meter® Professional Connectivity Solutions.

THERE'S STRENGTH IN OUR NUMBERS.

When you add up the numbers on our scale, it's not surprising that the 565KL-C is your choice for accuracy and reliability.

**5 MINUTE
WEIGHT HOLD**

**1 GRAM
RESOLUTION**

**120+ LOCAL
FIELD SPECIALISTS**

**#1 TRUSTED
BRAND**

**100 YEARS
OF SERVICE**

Health o meter®
Professional

Health o meter® Professional 565KL-C

weigh**easier**®



TECHNICAL SPECIFICATIONS

Capacity

44 lb / 20 kg

Resolution

0 – 44 lb x 0.1 oz / 0 – 6,000 g x 1 g /
6,001 g – 20,000 g x 2 g

EMR Connectivity

Via optional Pelstar®
wireless technology

Display

Digit Size: 3/4" / 19 mm
Screen Size: 3" x 1" LCD display

Power Source

100-240V Adapter (69-00001) included,
6 AA batteries not included

Functions

LB/ KG conversion
LB/ KG lock out
Zero
Tare
Pre-Tare
Hold / Release
Reweigh
Recall
Auto zero
Auto off
Motion detection
Auto-weight hold
Audible/Mute sound option
Variable auto-off time
Everlock® - permanent weight lock

Tray Size (w x d x h)

24 7/8" x 9 7/8" x 2 3/4"
(632 mm x 250 mm x 70 mm)

Product Footprint

26 1/4" x 13" x 5 1/4"
(667 mm x 330 mm x 134 mm)

Product Weight

9 lb / 4.1 kg

Measuring Tape

Range: 0 – 24 5/8" / 0 cm – 62.5 cm
Graduation: 1/16" / 1 mm

Portability

2210CART- optional rolling cart

Warranty

2-year limited manufacturer warranty

ORDERING INFORMATION

MODEL	DESCRIPTION
565KL-C	High resolution digital baby scale, LB & KG, Canada only
565KG-C	High resolution digital baby scale, KG only, Canada only
565EHR	Digital height rod for 565 series of scales
2210CART	Rolling scale/ medical device cart for 565 series of scales
SS-565KL	ScaleSurance® 2 year extended warranty for 565 series*
69-00001	Adapter included, item # for replacement or spare
BTPRINTERKIT-01	Wireless thermal printer kit

Scale also includes connectivity capabilities. To learn about connectivity options,
visit homscal.es.com/ConnectivitySolutions

Optional digital height rod
available



To learn more, call **1-800-815-6615**
or visit homscal.es.com/565KL-C

Health o meter®
Professional