

## THE 565KL

# HIGH-RESOLUTION DIGITAL BABY SCALE

The Health o meter® Professional 565KL neonatal/pediatric scale is designed to meet the demands of healthcare professionals in critical settings such as the NICU, nursery, delivery room, and pediatric practices. Its exceptional accuracy is especially crucial when caring for the tiniest patients.

The gently curved cradle provides comfort for infants, while the large angled display ensures caregivers can easily read weight measurements. Adjustable feet guarantee a level weighing surface, and the built-in measuring tape allows simultaneous weight and height assessment. Enhanced motion detection ensures quick locking even with an active infant. With a high 1-gram resolution, the scale accommodates a capacity of 20 kg (44 lb).

The 565KL's user-friendly digital interface offers a suite of features that enhance versatility in managing neonatal/pediatric patients. Notable functions include reweigh, hold/release, tare, recall, and a 5-minute weight hold. Additionally, the scale includes zero, pre-tare, KG/LB conversion, KG/LB/GM lockout, auto-zero, and auto-off features. The exclusive Everlock® function allows the measuring units (LB or KG) to be permanently locked on the scale.

The scale provides easy integration and seamless communication with today's EMR applications. The 565KL offers the option to incorporate wireless connectivity for interfacing with computers, EMR software, or other electronic devices via various Health o meter® Professional Connectivity Solutions.



## THERE'S STRENGTH IN OUR NUMBERS.

When you add up the numbers on our scale, it's not surprising that the 565KL is your choice for accuracy and reliability.

**5 MINUTE  
WEIGHT HOLD**

**1 GRAM  
RESOLUTION**

**120+ LOCAL  
FIELD SPECIALISTS**

**#1 TRUSTED  
BRAND**

**100 YEARS  
OF SERVICE**

**Health o meter®**  
Professional

Health o meter® Professional 565KL

weigh**easier**®

# TECHNICAL SPECIFICATIONS

## Capacity

44 lb / 20 kg

## Resolution

0 – 44 lb x 0.1 oz / 0 – 6,000 g x 1 g /  
6,001 g – 20,000 g x 2 g

## EMR Connectivity

Via optional Pelstar®  
wireless technology

## Display

Digit Size: 3/4" / 19 mm  
Screen Size: 3" x 1" LCD display

## Power Source

100-240V Adapter (69-00001) included,  
6 AA batteries not included

## Functions

LB/ KG conversion  
LB/ KG lock out  
Zero  
Tare  
Pre-Tare  
Hold/ Release  
Reweigh  
Recall  
Auto zero  
Auto off  
Motion detection  
Auto-weight hold  
Audible/Mute sound option  
Variable auto-off time  
Everlock® - permanent weight lock

## Tray Size (w x d x h)

24 7/8" x 9 7/8" x 2 3/4"  
(632 mm x 250 mm x 70 mm)

## Product Footprint

26 1/4" x 13" x 5 1/4"  
(667 mm x 330 mm x 134 mm)

## Product Weight

9 lb / 4.1 kg

## Measuring Tape

Range: 0 – 24 5/8" / 0 cm – 62.5 cm  
Graduation: 1/16" / 1 mm

## Portability

**2210CART**- optional rolling cart

## Warranty

2-year limited manufacturer warranty

# ORDERING INFORMATION

MODEL	DESCRIPTION
565KL	High resolution digital baby scale, LB & KG
565KG	High resolution digital baby scale, KG only
565EHR	Digital height rod for 565 series of scales
2210CART	Rolling scale/ medical device cart for 565 series of scales
SS-565KL	ScaleSurance® 2 year extended warranty for 565 series*
69-00001	Adapter included, item # for replacement or spare
BTPRINTERKIT-01	Wireless thermal printer kit

Scale also includes connectivity capabilities. To learn about connectivity options, visit [homscales.com/ConnectivitySolutions](https://homscales.com/ConnectivitySolutions)

Optional digital height rod available



To learn more, call **1-800-815-6615**  
or visit [homscales.com/565KL](https://homscales.com/565KL)

**Health o meter®**  
Professional

\*Not available in all countries.