

THE 2600ADAKIT

ADA COMPLIANCE KIT FOR 2600 SERIES OF SCALES

The Health o meter® Professional 2600ADAKIT is an optional ADA Compliance Kit for the 2600 series of scales. Following Americans with Disabilities Act (ADA) guidelines, the kit includes new ramps, new feet, and platform and ramp edge guards with anti-slip caution tape on the ramps to provide enhanced safety. The 2600ADAKIT is easy to install on any existing 2600 series scale, making it an ideal choice for doctor's offices, hospitals, clinics, and acute or long-term care facilities wanting to comply with the latest MDE standards.



2600ADAKIT for 2600 Series of Scales, Designed to Meet the US Access Boards Standards:

Medical facilities operated by state and local governments, also known as public entities, must comply with the MDE Standards for Accessible Medical Diagnostic Equipment, which is enforced by the Department of Justice.

- 32" wide x 40" deep platform
- 2" edge protection on platform
- 2" edge protection on ramps
- 1:8 running slope

ORDERING INFORMATION

MODEL	DESCRIPTION
2600ADAKIT	ADA Compliance Kit For 2600 Series Scales

Health o meter®
Professional

Health o meter® Professional 2600ADAKIT

TECHNICAL SPECIFICATIONS

Edge Guard Size (w x d x h)

40" x 1/8" x 2"
(1016 mm x 3 mm x 51 mm)

Ramp Size

32 1/4" x 16"
(819 mm x 406 mm)

Platform Size with 2600KL (w x d x h)

40" x 32 1/4" x 1 7/8"
(1016 mm x 819 mm x 48 mm)

Product Footprint with 2600KL

77 1/4" x 43 7/8" x 48 3/4"
(1962 mm x 1114 mm x 1238 mm)

Product Weight with 2600KL

126 lb / 57 kg

Warranty

2-year limited manufacturer warranty

2600ADAKIT Compatible Scales

2600KL
2600KG
2600KL-BT
2600KG-BT
2610KL
2610KG
2610KL-BT
2610KG-BT

weigh**easier**®

To learn more, call **1-800-815-6615**
or visit **homscales.com/2600ADAKIT**