


## THE 2620KL

# DIGITAL WHEELCHAIR SCALE WITH DUAL RAMPS AND EXTRA-LARGE PLATFORM, ADA COMPLIANT

The all-digital 2620KL is engineered with enhanced safety features to ensure a comfortable weighing experience for patients who use wheelchairs. Following Americans with Disabilities Act (ADA) compliance, the dual ramps adhere to a 1:8 slope ratio. The ramps and extra-large 32" x 40" platform feature a 2" edge protection to prevent roll-offs, and anti-slip caution tape on the ramps provide enhanced safety. The platform folds for storage, and built-in casters are included for added convenience and mobility. The scale has a weight capacity of up to 1,000 pounds (454 kilograms). Designed to meet U.S. Access Board Standards, the 2620KL is ideal for doctor's offices, hospitals, clinics, and acute or long-term care facilities. Additional features include an antimicrobial keypad, a swiveling display head with 180° rotation, and a large non-skid platform.

The scale features a high-contrast color screen with user-friendly functions like BMI, hold/release, reweigh, recall, zero, KG/LB, tare, and pre-tare. The Everlock® function permanently locks the measuring unit (LB or KG). With a USB interface, it seamlessly integrates with EMR systems and connects to devices via Health o meter® Professional Connectivity Solutions.



## THERE'S STRENGTH IN OUR NUMBERS.

When you add up the numbers on our scale, it's not surprising that the 2620KL provides a safe, accessible weighing experience.

**100%**  
ADA COMPLIANT

**2"** EDGE PROTECTION  
ON RAMPS  
AND PLATFORM

**32" X 40"** LOW  
PROFILE PLATFORM

**1:8** RUNNING  
SLOPE

**2620KL Wheelchair Scale, Designed to Meet the US Access Boards Standards:** Medical facilities operated by state and local governments, also known as public entities, must comply with the MDE Standards for Accessible Medical Diagnostic Equipment, which is enforced by the Department of Justice.

- 32" wide x 40" deep platform
- 2" edge protection on platform
- 2" edge protection on ramps
- 1:8 running slope

**Health o meter®**  
Professional

Health o meter® Professional 2620KL

**#1** TRUSTED  
BRAND

weigh**easier**®

# TECHNICAL SPECIFICATIONS

## Capacity

1000 lb / 454 kg

## Resolution

0.2 lb / 0.1 kg

## EMR connectivity

Via USB, optional Pelstar® wireless technology

## Display

Digit Size: 3/8" / 10 mm  
Screen Size: 2 3/4" x 1 1/2" High-Contrast  
Color TFT-LCD Screen

## Power Source

100-240V adapter (ADPT30) included, 6 D batteries not included

## Functions

LB/KG conversion  
LB/KG lock out  
BMI  
Zero  
Tare  
Pre-tare  
Hold/release  
Reweigh  
Recall  
Auto zero  
Auto off  
Time/date  
Variable auto off time  
Audible/mute sound option  
Parent/child weighing function  
Everlock®

## Platform Size (w x d x h)

40" x 32 1/4" x 1 7/8"  
(1016 mm x 819 mm x 48 mm)

## Ramp Size

32 1/4" x 16"  
(819 mm x 406 mm)

## Edge Guard Size (w x d x h)

40" x 1 1/8" x 2"  
(1016 mm x 3 mm x 51 mm)

## Product Footprint

77 1/4" x 43 7/8" x 48 3/4"  
(1962 mm x 1114 mm x 1238 mm)

## Product Weight

126 lb / 57 kg

## Portability

4 wheels (2 casters)  
for ease of movement

## Warranty

2-year limited manufacturer warranty

# ORDERING INFORMATION

MODEL	DESCRIPTION
2620KL	Digital Wheelchair Scale with Dual Ramps and Extra-Large Platform, ADA Compliant, LB & KG*
2620KG	Digital Wheelchair Scale with Dual Ramps and Extra-Large Platform, ADA Compliant, KG only*
SS-2620KL	ScaleSurance® 2 year extended warranty for 2620 series*
ADPT30	Adapter included, item # for replacement or spare
PLUGSET30	Plug adapters for use outside North America
BTPRINTERKIT-01	Wireless thermal printer kit

Scale also includes connectivity capabilities. To learn about connectivity options, visit [homscales.com/ConnectivitySolutions](https://homscales.com/ConnectivitySolutions)



Platform Folds for Mobility and Storage

To learn more, call **1-800-815-6615**  
or visit [homscales.com/2620KL](https://homscales.com/2620KL)



\*Not available in all countries.

**Health o meter®**  
Professional