


Product Specification Guide

Stand-On Scales

	499KL	500KL	502KL	599KL	600KL/ELEVATE
					
FEATURES					
Capacity	460 lb / 209 kg	660 lb / 300 kg	660 lb / 300 kg	660 lb / 300 kg	660 lb / 300 kg
Resolution	0.2 lb / 0.1 kg	0.2 lb / 0.1 kg	0.2 lb / 0.1 kg	0.2 lb / 0.1 kg	0.2 lb / 0.1 kg
Height Rod Included		•	•		•
Display Size	1" / 25 mm	1" / 25 mm	1" / 25 mm	1" / 25 mm	1" / 25 mm
Wheels	•	•	•	•	•
EMR Capable		Optional	Optional	•	•
Power Adapter Included				•	•
Batteries Included	•	•	•		
Warranty (years)	2	2	2	2	2
Wireless Connectivity		Optional	Optional	Optional	Optional
PRODUCT DIMS (w x d x h)					
Platform Size	13 ¾" x 16 ½" x 2 ¾"	13 ¾" x 16 ½" x 2 ¾"	13 ¾" x 16 ½" x 2 ¾"	14 ¼" x 14 ¼" x 2 ⅝"	14 ¼" x 14 ¼" x 2 ⅝"
Product Footprint	14 ⅞" x 21 ⅞" x 43 ½"	14 ⅞" x 21 ⅞" x 53 ⅝"	14 ⅞" x 21 ⅞" x 53 ⅝"	14 ¾" x 23 ¼" x 42 ½"	14 ¾" x 24 ⅞" x 54 ⅝"
Product Weight	21 lb / 9 kg	21 lb / 9 kg	22.4 lb / 10.2 kg	20 lb / 9 kg	28 lb / 13 kg
FUNCTIONS					
LB/KG Conversion	•	•	•	•	•
LB/KG Lockout	•	•	•	•	•
Hold/ Release	•	•	•	•	•
BMI	•	•	•	•	•
Zero	•	•	•	•	•
Tare	•	•	•	•	•
Pretare		•	•		
Auto Zero	•	•	•	•	•
Auto Off	•	•	•	•	•
Reweigh		•	•		
Recall		•	•		
Everlock®	•	•	•		
Parent/ Child Weighing	•	•	•	•	•
Variable Auto Off Time	•	•	•		
Audible/ Mute Sound	•	•	•		
Enabling Unit View		•	•		
ACCESSORIES (Optional)					
Power Adapter	ADPT31	ADPT31	ADPT31	ADPT50	ADPT50
Height Rod	499KLROD			499KLROD	
Wrap-Around Handlebar	499HB	500HB	500HB		
Wireless Kit		•	•	•	•
Wireless Printer		•	•	•	•

	1100KL	1110KL	2101KL	3001KL-AM	3001KL-AMX	3105KL-AM
						
FEATURES						
Capacity	1000 lb / 454 kg	1000 lb / 454 kg	1000 lb / 454 kg	1000 lb / 454 kg	1000 lb / 454 kg	1000 lb / 454 kg
Resolution	0.2 lb / 0.1 kg	0.2 lb/ 0.1 kg	0.2 lb / 0.1 kg	0.2 lb / 0.1 kg	0.2 lb / 0.1 kg	0.2 lb / 0.1 kg
Display Size	1 ½" / 38 mm	1 ½" / 38 mm	1 ½" / 38 mm	1 ½" / 38 mm	1 ½" / 38 mm	1 ½" / 38 mm
Wheels	•		•	•	•	•
EMR Capable	•	•	•	•	•	•
Power Adapter Included	•	•	•	•	•	
Batteries Included						•
Warranty (years)	2	2	2	2	2	2
Antimicrobial				•	•	•
Antimicrobial Keypad	•	•	•	•	•	•
Wireless Connectivity	Optional	Optional	Optional	Optional	Optional	Optional
PRODUCT DIMS (w x d x h)						
Platform Size	22" x 15 ¾" x 2 ½"	22" x 15 ¾" x 2 ⅝"	26" x 22" x 2 ¾"	22" x 15 ¾" x 2 ⅝"	22" x 15 ¾" x 2 ⅝"	19 ¾" x 21 ¼" x 2 ⅞"
Product Footprint	31 ¾" x 33 ½" x 49 ¾"	22" x 23 ½" x 4 ¾"	31" x 27" x 50"	22" x 34 ½" x 52"	22" x 34 ½" x 52"	23 ¼" x 28" x 52 ¼"
Product Weight	52 lb / 23 kg	27 lb / 12 kg	66 lb / 30 kg	42 lb / 19 kg	43 lb / 19.5 kg	86 lb / 39 kg
FUNCTIONS						
LB/KG Conversion	•	•	•	•	•	•
LB/KG Lockout	•	•	•	•	•	•
Hold/ Release	•	•	•	•	•	•
BMI	•	•	•	•	•	•
Zero	•	•	•	•	•	•
Tare	•	•	•	•	•	•
Pretare	•	•	•	•	•	•
Auto Zero	•	•	•	•	•	•
Auto Off	•	•	•	•	•	•
Reweigh	•	•	•	•	•	•
Recall	•	•	•	•	•	•
Everlock®	•	•	•	•	•	•
Parent/ Child Weighing	•	•	•	•	•	•
Time/ Date	•	•	•	•	•	•
Variable Auto Off Time	•	•	•	•	•	•
Audible/ Mute Sound	•	•	•	•	•	•
Enabling Unit View	•	•	•	•	•	•
ACCESSORIES (Optional)						
Power Adapter	ADPT30	ADPT30	ADPT30	ADPT30	ADPT30	ADPT30
Height Rod	PROPLUSROD or 245EHR-3105	201HR-1110 or 245EHR-1110	STROD or 245EHR-2101	3001HR or 245EHR-3001		201HR-3105 or 245EHR-3105
Wireless Kit	•	•	•	•	•	•
Wireless Printer	•	•	•	•	•	•

	349KLX	498KL	752KL	2842KL
FEATURES				
	Capacity	400 lb / 180 kg	660 lb / 300 kg	660 lb / 300 kg
	Resolution	0.2 lb / 0.1 kg	0.2 lb / 0.1 kg	0.2 lb / 0.1 kg
	Display Size	1" / 25 mm	1" / 25 mm	1" / 25 mm
	EMR Capable		Optional	•
	Power Adapter Included		•	•
	Batteries Included	•	•	
	Warranty (years)	2	2	2
	Wireless Connectivity		Optional	Optional
PRODUCT DIMS (w x d x h)				
Platform Size	12 ½" x 12" x 1 ⅞"	13 ¾" x 16 ½" x 2 ⅝"	14 ¼" x 14 ¼" x 2 ⅝"	22 ¼" x 42" x 2 ½"
Product Footprint	12 ½" x 12" x 1 ⅞"	13 ¾" x 16 ½" x 2 ⅝"	14 ½" x 14 ½" x 2 ¾"	22 ¼" x 42" x 2 ½"
Product Weight	10 lb / 4.5 kg	14 lb / 6.4 kg	14 lb / 6.4 kg	53 lb / 24 kg
FUNCTIONS				
LB/KG Conversion	•	•	•	•
LB/KG Lockout	•	•	•	
Hold/ Release	•	•	•	
BMI		•	•	
Zero	•	•	•	•
Tare	•	•	•	•
Pretare		•		
Auto Zero	•	•	•	
Auto Off	•	•	•	
Reweigh	•	•		
Recall				
Everlock®	•	•		
Parent/ Child Weighing	•	•	•	
Variable Auto Off Time	•	•		
Audible/ Mute Sound		•		
Enabling Unit View	•	•		
ACCESSORIES (Optional)				
Power Adapter	ADPT40	ADPT31	ADPT50	ADPT31
Carrying Case	64771	64771	64771	
Wireless Kit		•	•	
Wireless Printer		•	•	

Mechanical Stand-On Scales

	400KL	402KL	402LB	402KG SPECIAL ORDER ONLY	450KL
FEATURES					
	Capacity	390 lb / 180 kg	390 lb / 180 kg	400 lb	180 kg
	Graduation	¼ lb / 100 g	¼ lb / 100 g	¼ lb / 100 g	¼ lb / 100 g
	Height Rod Included		•	•	•
	Rotating Poise Bars				•
	Engraved Numerals				•
	Warranty (years)	10	10	10	10
PRODUCT DIMS (w x d x h)					
Platform Size	10 ½" x 14" x 3 ¼"	10 ½" x 14" x 3 ¼"	10 ½" x 14" x 3 ¼"	10 ½" x 14" x 3 ¼"	10 ½" x 14" x 3 ¼"
Product Footprint	20 ¼" x 10 ⅝" x 57 ½"	20 ¼" x 10 ⅝" x 57 ½"	20 ¼" x 10 ⅝" x 57 ½"	20 ¼" x 14" x 57 ½"	20 ¼" x 14" x 57 ½"
Product Weight	24 lb / 11 kg	28 lb / 13 kg	28 lb / 13 kg	28 lb / 13 kg	31 lb / 14 kg
ACCESSORIES (Optional)					
Wheels	55000	55000	55000	55000	55000
Height Rod	402KLROD				
Counterweights	55090	55090			Included

Pediatric Scales



	522KL	524KL	549KL	550KL	553KL
					
FEATURES					
Capacity	50 lb / 23 kg	50 lb / 23 kg	44 lb / 20 kg	44 lb / 20 kg	44 lb / 20 kg
Resolution	0-20 lb / 0.2 oz; 20-50 lb / 0.5 oz; 0-9 kg / 5 g; 9-23 kg / 10 g	0-20 lb / 0.2 oz; 20-50 lb / 0.5 oz; 0-9 kg / 5 g; 9-23 kg / 10 g	0.5 oz / 10 g	0.5 oz / 10 g	0.5 oz / 10 g
Measuring Tape Included	•				•
Display Size	1" / 25 mm	1" / 25 mm	5⁄8" / 15 mm	½" / 13 mm	1" / 25 mm
EMR Capable	Optional	Optional			Optional
Batteries Included	•	•	•	•	•
Warranty (years)	2	2	2	2	2
Wireless Connectivity	Optional	Optional			Optional
PRODUCT DIMS (w x d x h)					
Tray/ Seat Size	24 ⅙" x 14 ½" x 2 ⅝"	24 ¼" x 14 ¼" x 18 ½"	19 ⅜" x 12 ⅙" x 2"	23 ½" x 8" x 2"	24 ⅞" x 13" x 1 ⅝"
Product Footprint	24 ⅙" x 21" x 23 ⅙"	24 ⅙" x 21" x 23 ⅙"	19 ⅜" x 13 ½" x 4 ⅝"	23 ½" x 13 ¼" x 3 ⅝"	24 ⅞" x 16 ½" x 3 ⅞"
Product Weight	20 lb / 9 kg	21 lb / 10 kg	4 lb / 2 kg	8 lb / 4 kg	10 lb / 5 kg
FUNCTIONS					
LB/KG Conversion	•	•	•		•
LB/KG Lockout	•	•	•		•
Hold/ Release	•	•			•
Zero	•	•			•
Tare	•	•	•	•	•
Pretare	•	•			•
Auto Zero	•	•	•		•
Auto Off	•	•	•	•	•
Reweigh	•	•			•
Recall	•	•			•
Everlock®	•	•	•		•
Variable Auto Off Time	•	•	•		•
Audible/ Mute Sound	•	•			
Enabling Unit View	•	•	•		•
Motion Detection	•	•	•		•
Auto Weight Hold	•	•	•	•	•
ACCESSORIES (Optional)					
Power Adapter	ADPT31	ADPT31			ADPT31
Carrying Case			64769	553CASE	553CASE
Height Rod	522EHR or 522HR			eBABYHR + 550EHRBRACKET	553HR or 553EHR
Cart			2210CART	2210CART	2210CART
Wireless Kit	•	•			•
Wireless Printer	•	•			•

	565KL	2210KL-AM	2210KL4-AM	8320KL
				
FEATURES				
Capacity	44 lb / 20 kg	45 lb / 20 kg / 20,000 gm	45 lb / 20 kg / 20,000 gm	33 lb / 15 kg
Resolution	0 – 44 lb x 0.1 oz / 0 – 6,000 g x 1 g / 6,001 g – 20,000 g x 2 g	1 g (to 10,000 g); 5 g (10,001 to 20,000 g)	1 g (to 10,000 g); 5 g (10,001 to 20,000 g)	0.5 oz/ 10 g
Measuring Tape Included	•	•	•	
Display Size	¾" / 19 mm	1 ⅞" / 28 mm	1 ⅞" / 28 mm	¾" / 19 mm
EMR Capable	Optional	•	•	
Batteries Included	•	•	•	•
Warranty (years)	2	2	2	2
Antimicrobial		•	•	
Wireless Connectivity	Optional	Optional	Optional	
PRODUCT DIMS (w x d x h)				
Tray/ Seat Size	24 ⅞" x 9 ⅞" x 2 ¾"	24" x 15 ¼" x 2 ¾"	26 ½" x 15 ¼" x 2 ¾"	23 ½" x 11" x ⅞"
Product Footprint	26 ¼" x 13" x 5 ¼"	24 ¾" x 17 ¾" x 6"	26 ½" x 17 ¾" x 6"	29 ⅞" x 12 ⅝" x 2 ⅞"
Product Weight	9 lb / 4.1 kg	28 lb / 13 kg	28 lb / 13 kg	5 lb / 3 kg
FUNCTIONS				
LB/KG Conversion	•	•	•	•
LB/KG Lockout	•	•	•	
Hold/ Release	•	•	•	
Zero	•	•	•	•
Tare	•	•	•	
Pretare	•	•	•	
Auto Zero	•	•	•	•
Auto Off	•	•	•	•
Reweigh	•	•	•	
Recall	•	•	•	
Everlock®	•	•	•	
Variable Auto Off Time	•	•	•	
Audible/ Mute Sound	•	•	•	
Enabling Unit View	•	•	•	
Motion Detection	•	•	•	
Auto Weight Hold	•	•	•	
Breast Milk Intake		•	•	
ACCESSORIES (Optional)				
Power Adapter	69-00001	ADPT30	ADPT30	
Carrying Case				64771
Height Rod	565EHR	2210EHR		
Cart	2210CART	2210CART-AC	2210CART-AC	
Wireless Kit	•	•	•	
Wireless Printer	•	•	•	

	1522KL	1524KL	386S-01	386KGS-01	BABYSLINGSCALE
					
FEATURES					
Capacity	135 lb / 62 kg	130 lb / 65 kg	50 lb	25 kg	25 kg
Graduation	1 oz / 20 g	1 oz / 20 g	¼ lb	100 g	100 g
Measuring Tape Included	•		•	•	
Dial Size			3" / 76 mm	3" / 76 mm	5 ½" / 140 mm
Warranty (years)	2	2	1	1	1
PRODUCT DIMS (w x d x h)					
Tray/ Seat Size	24" x 16 ½" x 2 ⅞"	24 ½" x 14 ¼" x 18 ½"	19 ⅝" x 12 ⅝" x 3 ¾"	19 ⅝" x 12 ⅝" x 3 ¾"	
Product Footprint	24" x 20" x 22 ½"	24 ½" x 21" x 22 ½"	19 ⅝" x 14 ½" x 5 ⅛"	19 ⅝" x 14 ½" x 5 ⅛"	12" x 43" x 1 ⅞"
Product Weight	25 lb / 11 kg	32 lb / 15 kg	10 lb / 5 kg	10 lb / 5 kg	2 lb / 1 kg

Wheelchair & Chair Scales

	2500KL	2600KL	2610KL	2620KL
				
FEATURES				
Capacity	1000 lb / 454 kg	1000 lb / 454 kg	1000 lb / 454 kg	1000 lb / 454 kg
Resolution	0.2 lb / 0.1 kg	0.2 lb / 0.1 kg	0.2 lb / 0.1 kg	0.2 lb / 0.1 kg
Display Size	1 ½" / 38 mm	1 ½" / 38 mm	1 ½" / 38 mm	1 ½" / 38 mm
Wheels	•	•	•	•
EMR Capable	•	•	•	•
Power Adapter Included	•	•	•	•
Warranty (years)	2	2	2	2
Antimicrobial				
Antimicrobial Keypad	•	•	•	•
Wireless Connectivity	Optional	Optional	Optional	Optional
ADA Compliant	•			
PRODUCT DIMS				
Platform Size	29 ⅝" x 28 ¼" x 2 ⅞"	36" x 32 ¼" x 2 ½"	36" x 32 ¼" x 2 ½"	40" x 32 ¼" x 1 ⅞"
Product Footprint	40 ¼" x 34 ½" x 45 ½"	48 ⅛" x 43 ⅞" x 49 ⅜"	56 ¾" x 43 ⅞" x 49 ⅜"	77 ¼" x 43 ⅞" x 48 ¾"
Product Weight	81 lb / 37 kg	93 lb / 42 kg	98 lb / 44 kg	126 lb / 57 kg
FUNCTIONS				
LB/ KG Conversion	•	•	•	•
LB/ KG Lockout	•	•	•	•
Hold/ Release	•	•	•	•
BMI	•	•	•	•
Zero	•	•	•	•
Tare	•	•	•	•
Pretare	•	•	•	•
Auto Zero	•	•	•	•
Auto Off	•	•	•	•
Reweigh	•	•	•	•
Recall	•	•	•	•
Everlock®	•	•	•	•
Time/ Date	•	•	•	•
Parent/ Child Weighing	•	•	•	•
Variable Auto Off Time	•	•	•	•
Audible/ Mute Sound	•	•	•	•
Enabling Unit View	•	•	•	•
ACCESSORIES (Optional)				
Power Adapter	ADPT30	ADPT30	ADPT30	ADPT30
Wireless Kit	•	•	•	•
Wireless Printer	•	•	•	•
Ramp	B2500RAMP	B2600RAMP		
ADA Compliance Upgrade Kit			2600ADAKIT	2600ADAKIT

	2400KL	2650KL	2700KL	2900KL-AM
				
FEATURES				
Capacity	800 lb / 363 kg	1000 lb / 454 kg	1000 lb / 454 kg	1000 lb / 454 kg
Resolution	0.2 lb / 0.1 kg	0.2 lb / 0.1 kg	0.2 lb / 0.1 kg	0.2 lb / 0.1 kg
Display Size	1 ½"/ 38 mm	1 ½"/ 38 mm	1 ½"/ 38 mm	1 ½"/ 38 mm
Wheels	•	•	•	
EMR Capable	•	•	•	•
Power Adapter Included	•	•	•	•
Warranty (years)	2	2	2	2
Antimicrobial				•
Antimicrobial Keypad	•	•	•	•
Wireless Connectivity	Optional	Optional	Optional	Optional
PRODUCT DIMS (w x d x h)				
Platform Size	25 ¾" - 50 ½"	31 ½" x 31 ½" x 2"	43" x 42" x 2 ½"	29 ½" x 32 ½" x 2 ¾"
Product Footprint	51 ¼" x 40 ½" x 3 ½"	51 ¾" x 33 ¾" x 5 ¼"	69 ⅞" x 51 ½" x 3 ⅝"	58" x 37" x 37"
Product Weight	40 lb / 18 kg	82 lb / 37 kg	187 lb / 85 kg	91 lb / 41 kg
FUNCTIONS				
LB/ KG Conversion	•	•	•	•
LB/ KG Lockout	•	•	•	•
Hold/ Release	•	•	•	•
BMI	•	•	•	•
Zero	•	•	•	•
Tare	•	•	•	•
Pretare	•	•	•	•
Auto Zero	•	•	•	•
Auto Off	•	•	•	•
Reweigh	•	•	•	•
Recall	•	•	•	•
Everlock®	•	•	•	•
Time/ Date	•	•	•	•
Parent/ Child Weighing	•	•	•	•
Variable Auto Off Time	•	•	•	•
Audible/ Mute Sound	•	•	•	•
Enabling Unit View	•	•	•	•
ACCESSORIES (Optional)				
Power Adapter	ADPT30	ADPT30	ADPT30	ADPT30
Wireless Kit	•	•	•	•
Wireless Printer	•	•	•	•

	2500CKL	2595KL	594KL	445KL
				
FEATURES				
Capacity	1000 lb / 454 kg	600 lb / 270 kg	600 lb / 270 kg	440 lb / 200 kg
Resolution	0.2 lb / 0.1 kg	0.2 lb / 0.1 kg	0.2 lb / 0.1 kg	0.2 lb / 0.1 kg
Display Size	1 ½"/ 38 mm	1 ⅝"/ 30 mm	1 "/ 25 mm	
Wheels	•	•	•	•
EMR Capable	•			
Power Adapter Included	•	•	•	
Warranty (years)	2	2	2	2
Antimicrobial Keypad	•			
Wireless Connectivity	Optional	Optional		
PRODUCT DIMS (w x d x h)				
Platform Size	29 ¾" x 28 ¼" x 2 ¾"			
Seat Size	18" x 15"	18 ¼" x 15 ¾"	18" x 18 ½"	18" x 17 ¼"
Product Footprint	40 ¼" x 34 ½" x 45 ½"	23 ½" x 37 ½" x 39 ⅞"	23 ¼" x 40 ⅞" x 37 ⅞"	24 ¾" x 31 ⅝" x 39 ⅞"
Product Weight	121 lb / 55 kg	84 lb / 38 kg	39 lb / 18 kg	60 lb / 27 kg
FUNCTIONS				
LB/ KG Conversion	•	•	•	
LB/ KG Lockout	•	•	•	
Hold/ Release	•	•	•	
BMI	•		•	
Zero	•	•	•	
Tare	•		•	
Pretare	•	•		
Auto Zero	•		•	
Auto Off	•	•	•	
Reweigh	•	•		
Recall	•	•		
Everlock®	•	•		
Time/ Date	•			
Parent/ Child Weighing	•			
Variable Auto Off Time	•	•		
Audible/ Mute Sound	•			
Enabling Unit View	•	•		
ACCESSORIES (Optional)				
Power Adapter	ADPT30	ADPT30	ADPT30	
Wireless Kit	•	•		
Wireless Printer	•			
Ramp	B2500RAMP			

Body Composition Scales



FEATURES

Capacity	1000 lb / 454 kg
Resolution	0.2 lb / 0.1 kg
Body Fat Accuracy	+/- 2%
Measurement Time	15 seconds
Power Adapter Included	•
6 Window Display	•
Limited Test Restrictions	•
OHMS Modulation Logic	•
Warranty (years)	1

MEASUREMENTS

Weight	•
% Body Fat	•
% Body Water	•
% Muscle Mass	•
BMI (OR BMI %)	•
REE/ BMR	•
Activity-Based Caloric Burn Rate	•

IMPACT ON ACCURACY/ PRE-TEST RESTRICTIONS

Recent Food/ Drink	OK to Test
Recent Exercise	OK to Test
Full Night Sleep	OK to Test
Recent Alcohol Intake	OK to Test
Full Bladder	OK to Test

REPORTS

Plain Paper	•
Illustrated	•
Off-the-Shelf Printer	•

PRODUCT DIMS (w x d x h)

Platform Size	14" x 16" x 2 7⁄8"
Product Footprint	18" x 23 ½" x 41 ¼"
Product Weight	27 lb / 12 kg

USAGE

	BCS-G6-adult	BCS-G6-limbs	BCS-G6-DUO	BCS-G61-upper
Adult Patients	•	•	•	•
Pediatric Patients			•	
Segmented (Incl. Arms & Legs)		•		
Upper Body Only				•
Shoes OK				•

Professional Home Care Scales

800KL

810KL

880KL

844KL

822KL



FEATURES

Capacity	397 lb / 180 kg	440 lb / 200 kg	500 lb / 220 kg	440 lb / 200 kg	400 lb / 180 kg
Resolution	0.2 lb / 0.1 kg	0.2 lb / 0.1 kg	0.2 lb / 0.1 kg	0.1 lb / 0.05 kg	0.1 lb / 0.05 kg
Display Size	1 ½" / 38 mm	1 ⅞" / 48 mm	1 ⅞" / 48 mm	1 ¼" / 32 mm	1 ½" / 38 mm
Batteries Included	•	•	•	•	•
Case Pack Quantity	4	1	4	2	2
Warranty (years)	1	1	1	1	1

PRODUCT DIMS (w x d x h)

Platform Size	12 ¼" x 12 ¼" x 1 ¾"	13" x 13" x 7⁄8"	17" x 13" x 1 ¼"	12 ⅝" x 12 ⅝" x 1 7⁄8"	14 ¼" x 11 ¾" x 1 ½"
Product Footprint	12 ¼" x 12 ¼" x 1 ¾"	13" x 13" x 7⁄8"	17" x 13" x 1 ¼"	12 ⅝" x 12 ⅝" x 1 7⁄8"	14 ¼" x 11 ¾" x 1 ¾"
Product Weight	4 lb / 2 kg	4 lb / 2 kg	8 lb / 4 kg	4 lb / 2 kg	4 lb / 2 kg

FUNCTIONS

LB/ KG Conversion		•			
Auto Zero	•	•	•	•	•
Auto Off	•	•	•	•	•
5 Second Weight Hold		•			

ACCESSORIES (Optional)

Carrying Case	64771		64769	64771	64771
---------------	-------	--	-------	-------	-------

	142KL	160KL	160LB	160KG
FEATURES				
	Capacity	330 lb / 150 kg	400 lb / 180 kg	400 lb
	Resolution	1 lb / 0.5 kg	1 lb / 1 kg	1 lb
	Dial Size	6 ½" / 165 mm	8" / 203 mm	8" / 203 mm
	Case Pack Quantity	2	2	2
	Warranty (years)	1	1	1
PRODUCT DIMS (w x d x h)				
Platform Size	11 ½" x 10 ¾" x 2 ¼"	11" x 12 ½" x 3"	11" x 12 ½" x 3"	11" x 12 ½" x 3"
Product Footprint	11 ¾" x 17 ⅝" x 3 ½"	11 ¾" x 18 ¼" x 7 ¼"	11 ¾" x 18 ¼" x 7 ¼"	11 ¾" x 18 ¼" x 7 ¼"
Product Weight	6 lb / 3 kg	9 lb / 4 kg	9 lb / 4 kg	9 lb / 4 kg


	100LB	100KG	175LB
FEATURES			
	Capacity	270 lb	120 kg
	Resolution	1 lb	1 lb
	Dial Size	2" / 51 mm	2" / 51 mm
	Case Pack Quantity	3	3
	Warranty (years)	1	1
PRODUCT DIMS (w x d x h)			
Platform Size	10 ¼" x 10" x 1 ⅞"	10 ¼" x 10" x 1 ⅞"	10 ½" x 9 ⅝" x 2"
Product Footprint	10 ¼" x 10" x 2"	10 ¼" x 10" x 2"	11 ⅝" x 9 ¾" x 2 ¼"
Product Weight	4 lb / 2 kg	4 lb / 2 kg	4 lb / 2 kg
ACCESSORIES (Optional)			
Carrying Case	64771	64771	64771

Specialty Scales

	PTS-1000KL	2000KL	3400KL	3401KL	YG500R	
FEATURES						
	Capacity	550 lb / 250 kg	400 lb / 181 kg	176 oz / 5000 g	176 oz / 5000 g	500 g
	Resolution/ Graduation	1 lb / 0.5 kg	0.2 lb / 0.1 kg	0.05 oz / 1 g	0.05 oz / 1 g	2 g
	Display/ Dial Size	1" / 25 mm	1 ½" / 38 mm	⅝" / 16 mm	⅝" / 16 mm	6" / 152 mm
	Wheels					•
	EMR Capable					•
	Power Adapater Included	•	•	•	•	
	Warranty (years)	2	2	2	2	1
Wireless Connectivity	Optional					
PRODUCT DIMS (w x d x h)						
Platform Size	71" x 21"	70 ½" x 30"	9 ½" x 11 ¾" x 2 ½"	7 ¾" x 9 ¾" x ⅝"	6 ½" x 4 ⅔"	
Product Footprint	71" x 27 ½" x 1"	70 ½" x 46 ½" x 60 ½"	10 ⅞" x 13" x 8"	10" x 11 ¾" x 5 ⅝"	6 ½" x 6 ⅞" x 8 ⅝"	
Product Weight	25 lb / 11.4 kg	104 lb / 47 kg	11 lb / 5 kg	10 lb / 4.5 kg	3 lb / 1 kg	
FUNCTIONS						
LB/ KG Conversion	•	•	•	•		
LB/ KG Lockout	•					
Hold/ Release	•	•				
BMI	•					
Zero		•	•	•		
Tare	•					
Pretare	•					
Auto Zero	•					
Auto Off	•	•	•	•		
Reweigh	•					
Recall	•					
Everlock®	•					
Time/ Date	•					
Variable Auto Off Time	•					
Audible/ Mute Sound	•					
Enabling Unit View	•					
ACCESSORIES (Optional)						
Power Adapter	69-00012	APT30	3400ADPT	3400ADPT		
Wireless Kit	•					
Wireless Printer	•					
Other			3400PAN (3 pack pans) 3400TRAY (3 pack trays)	3400PAN (3 pack pans) 3400TRAY (3 pack trays)		

Height Measuring Instruments

	201HR	205HR	PORTROD	WMROD
				
FEATURES				
Measurement Range	2" - 90" / 6 - 230 cm	0 - 90" / 0 - 230 cm	24" - 83" / 60 - 213 cm	11 7⁄8" - 84" / 30 - 213 cm
Graduation	1/16" / 1 mm	1/16" / 1 mm	1/16" / 1 mm	1/8" / 1 mm
Warranty (years)	2	2	2	2
PRODUCT DIMS (w x d x h)				
Product Footprint	2 ¼" x 11" x 54"	2 ⅝" x 8" x 60"	2 ¾" x 8 ¼" x 60"	1 ½" x 13 ½" x 80 ¼"
Product Weight	2 lb / 1 kg	2 lb / 1 kg	1 lb / .5 kg	4.2 lb / 2 kg

	221HR	280HR	BABYHR	EBABYHR
				
FEATURES				
Measurement Range	7" - 82.6" / 17.7 - 210 cm	3.9 - 43.3" / 10 - 110 cm	3" - 39" / 7 - 100 cm	13 ¾" - 31 ½" / 35 - 80 cm
Graduation	1/16" / 1 mm	1/4" / 5 mm	1/16" / 1 mm	1/16" / 1 mm
Warranty (years)	2	2	2	2
PRODUCT DIMS (w x d x h)				
Product Footprint	13 ¾" x 23" x 89"	55" x 12"	1 ½" x 11 ¼" x 41"	2 ¾" x 11 ½" x 24 ¾"
Product Weight	7.9 lb / 3.6 kg	.88 lb / .4 kg	1 lb / .5 kg	1 lb / .5 kg

Height Measuring Accessories for Platform Scales

	Height Rods for 1100KL Series	Height Rods for 1110KL Series	Height Rods for 2101KL Series	Height Rods for 3001KL-AM Series	Height Rods for 3105KL-AM Series
	PROPLUSROD	201HR-1110	STROD	3001HR	201HR-3105
					
FEATURES					
Mechanical Measurement Range	2" - 90" / 6 - 230 cm	2" - 90" / 6 - 230 cm	2" - 90" / 6 - 230 cm	2" - 90" / 6 - 230 cm	2" - 90" / 6 - 230 cm
Mechanical Measurement Graduation	1/16" / 1 mm	1/16" / 1 mm	1/16" / 1 mm	1/16" / 1 mm	1/16" / 1 mm
Warranty (years)	2	2	2	2	2
PRODUCT DIMS (w x d x h)					
Product Footprint	2 ¼" x ¾" x 54"	2 ¼" x ¾" x 54"	2 ¼" x ¾" x 54"	2 ¼" x ¾" x 54"	2 ¼" x ¾" x 54"
Product Weight	2 lb / 1 kg	2 lb / 1 kg	2 lb / 1 kg	2 lb / 1 kg	2 lb / 1 kg
	245EHR-1100	245EHR-1110	245EHR-2101	245EHR-3001	245EHR-3105
					
FEATURES					
Mechanical Measurement Range	11 ¾" - 79 ½" / 30 - 202 cm	11 ¾" - 79 ½" / 30 - 202 cm	11 ¾" - 79 ½" / 30 - 202 cm	11 ¾" - 79 ½" / 30 - 202 cm	11 ¾" - 79 ½" / 30 - 202 cm
Mechanical Measurement Graduation	1/16" / 1 mm	1/16" / 1 mm	1/16" / 1 mm	1/16" / 1 mm	1/16" / 1 mm
Digital Measurement Range	3 ft. 11.2" - 6 ft. 7.5" / 120 cm- 202 cm	3 ft. 11.2" - 6 ft. 7.5" / 120 cm- 202 cm	3 ft. 11.2" - 6 ft. 7.5" / 120 cm- 202 cm	3 ft. 11.2" - 6 ft. 7.5" / 120 cm- 202 cm	3 ft. 11.2" - 6 ft. 7.5" / 120 cm- 202 cm
Digital Measurement Resolution	0.1" / 1 mm	0.1" / 1 mm	0.1" / 1 mm	0.1" / 1 mm	0.1" / 1 mm
Warranty (years)	2	2	2	2	2
PRODUCT DIMS (w x d x h)					
Product Footprint	2 ¼" x 2" x 48"	2 ¼" x 2" x 48"	2 ¼" x 2" x 48"	2 ¼" x 2" x 48"	2 ¼" x 2" x 48"
Product Weight	4 lb / 2 kg	4 lb / 2 kg	4 lb / 2 kg	4 lb / 2 kg	4 lb / 2 kg

	Height Rod for 400KL Series	Height Rod for 499KL & 599KL Series
	402KLROD	499KLROD
		
FEATURES		
Mechanical Measurement Range	11 7⁄8" - 84" / 30 - 213 cm	11 7⁄8" - 84" / 30 - 213 cm
Mechanical Measurement Graduation	1/8" / 1 mm	1/8" / 1 mm
Warranty (years)	2	2
PRODUCT DIMS (w x d x h)		
Product Footprint	1 ½" x 13 ½" x 80 ¼"	1 ½" x 13 ½" x 80 ¼"
Product Weight	4.2 lb / 2 kg	4.2 lb / 2 kg

Height Measuring Accessories for Pediatric Scales

	Height Rods for 522KL Series	Height Rods for 553KL Series	Height Rod for 550KL	Height Rod for 565KL Series	Height Rod for 2210KL-AM Series
	522EHR	553EHR	eBABYHR + 550EHRBRACKET	565EHR	2210EHR
					
FEATURES					
Digital Measurement Range	13 3⁄4" – 31 1⁄2" / 35 cm – 80 cm	13 3⁄4" – 31 1⁄2" / 35 cm – 80 cm	13 3⁄4" – 31 1⁄2" / 35 cm – 80 cm	13 3⁄4" – 31 1⁄2" / 35 cm – 80 cm	13 3⁄4" – 31 1⁄2" / 35 cm – 80 cm
Digital Measurement Resolution	1/16" / 0.1 cm	1/16" / 0.1 cm	1/16" / 0.1 cm	1/16" / 0.1 cm	1/16" / 0.1 cm
Warranty (years)	2	2	2	2	2
PRODUCT DIMS (w x d x h)					
Product Footprint	2 3⁄4" x 11 1⁄2" x 24 3⁄4"	2 3⁄4" x 11 1⁄2" x 24 3⁄4"	2 3⁄4" x 11 1⁄2" x 24 3⁄4"	2 3⁄4" x 11 1⁄2" x 24 3⁄4"	2 3⁄4" x 11 1⁄2" x 24 3⁄4"
Product Weight	2 lb / 1 kg	2 lb / 1 kg	2 lb / 1 kg	2 lb / 1 kg	2 lb / 1 kg
	522HR	553HR			
					
FEATURES					
Mechanical Measurement Range	3" - 39" / 7 - 100 cm	3" - 39" / 7 - 100 cm			
Mechanical Measurement Graduation	1/16" / 1 mm	1/16" / 1 mm			
Warranty (years)	2	2			
PRODUCT DIMS (w x d x h)					
Product Footprint	1 1⁄2" x 11 1⁄4" x 41"	1 1⁄2" x 11 1⁄4" x 41"			
Product Weight	1 lb / .5 kg	1 lb / .5 kg			